

New features in Mascot Server 2.7

Richard Jacob

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



New features in Mascot Server 2.7

- **Crosslinking, Quantitation Summaries and Modifications**
 - Are covered in full in separate talks
- **Protein FDR**
- **Fragment charge states**
- **Better filtering of modifications with complementary deltas in an error tolerant searches**
- **Combined search results in Mascot Daemon**
- **MGF exporting spectra in the original order**
- **Handle non-standard modification names in a spectral library**

The big new feature support for crosslinking. This and the Quantitation Summaries features and improvements to how we search for Modifications are covered in separate talks.

In this talk I will cover the following topics:

Crosslinking

- **Naturally occurring**
 - E.g. disulfide bond
- **Cleavable crosslinks**
 - Cleaved prior to analysis, e.g. irradiation or CID
 - Leaves behind a fragment at linked site
- **Intact crosslinks**
 - Survives MS/MS analysis
 - Spectrum has fragments from both peptides

We've added full support for a comprehensive range of crosslinking scenarios. This includes editors for the crosslinking agents and methods, fragmentation patterns in peptide_view and exports to xiVIEW for further visualization.

Crosslinking

88	6.1e-09	▶1	U K.SLQSNLTLEEIQERDAKLRK.E + K17<-Xlink:DSS[...]->K20
106	1.2e-10	▶1	U K.SLQSNLTLEEIQERDAKLRK.E + Xlink:DSS[W]
123	3e-12	▶1	U K.SLQSNLTLEEIQERDAKLRK.E + Xlink:DSS[W]
49	7.5e-05	▶1	U K.EMEEKLVKLRREGITLVRPEDK.K
59	3.2e-06	▶1	U K.TAVQKALDSLADAGK.I K5<-Xlink:DSS->K2 K.QKIYIAR.Q
34	0.00062	▶1	U K.TAVQKALDSLADAGK.I K5<-Xlink:DSS->K2 K.QKIYIAR.Q
46	5.4e-05	▶1*	[2:MND1_ARATH]LESDLQGSNK K10<-Xlink:DSS->K8 K.LQEQLQEKK.K

Looplinked peptide

Monolinks

Intralinks in HOP2_ARATH

Interlink to MND1_ARATH

MASCOT

: New features in Mascot 2.7

© 2020 Matrix Science

MATRIX SCIENCE

Mascot can report all the different types of links from looped, Inter and Intra links and incomplete mono links.

Modifications

Mascot Server 2.6

- 2 methods
- Less than 16 possible permutations:
- More than 16 possible permutations:
 - Tends to cluster modifications on adjacent modifiable sites
 - Often stops before 16 different permutations tested

Mascot Server 2.7

- Single, consistent, permutation method
- Controlled by 3 user definable parameters:
 - MaxPepNumVarMods
 - MaxPepNumModifiedSites
 - MaxPepModArrangements

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science




Mascot Server 2.7 has better site localisation by introducing separate caps on the number of distinct modifications in a PSM, the number of modified sites in a PSM, and the total number of arrangements to be tested.

Mascot Server 2.6 worked well but had limitations. If a peptide has less than 16 possible variable modification permutations, then there are no issues as all possibilities are tested. It's in the cases where the peptide has more than 16 possible permutations that issues can arise. The sliding window method tends to cluster modifications on adjacent modifiable sites, and it will often stop before 16 different permutations have been tested.

In Mascot 2.7 we've taken a different approach. We use a single, consistent permutation method – there's no switching between different methods. The new permutation iterator samples arrangements using a uniformly random scheme. The operation of this is controlled by 3 user definable settings.

By adjusting the parameters you can choose to reduce the search depth and speed up searches of lightly modified samples or increase the search depth where site localization matters or for highly modified samples.

Quantitation Summaries



MASCOT : New features in Mascot 2.7

© 2020 Matrix Science



The Quantitation Summaries feature in Mascot Daemon allows you to select one or more tasks and create a sample map.

Sample map

<input checked="" type="checkbox"/>	Task	Task name	Raw file path	Raw file name	Parameter set	Fastq(s)	Submitted	Fraction	Intensity LFG
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_1_3ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 17.2	1	3ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_2_3ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 18.0	2	3ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_3_3ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 18.4	3	3ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_4_7.5ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 19.2	1	7.5ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_5_7.5ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 20.0	2	7.5ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_6_7.5ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 20.4	3	7.5ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_7_10ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 21.1	1	10ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_8_10ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 22.0	2	10ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_8_10ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 22.4	3	10ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_10_15ng_270914.raw	weizmann.par	ecoli_protome...	06/08/2019 23.2	1	15ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_11_15ng_270914.raw	weizmann.par	ecoli_protome...	07/08/2019 00.0	2	15ng
<input checked="" type="checkbox"/>	524	Weizmann LFG	\kippy\matrasci	GEF1_Spiketh_230914_12_15ng_270914.raw	weizmann.par	ecoli_protome...	07/08/2019 00.4	3	15ng

Contaminant DB: [None] | Average [MC] | Settings | Save sample map | Save quantitation summary | Close

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



The Sample Map is used to annotate the list of result files with recognisable sample identifiers

MASCOT : New features in Mascot 2.7 © 2020 Matrix Science

And Daemon creates and saves the Quantitation Summary, a tab delimited text file which is the input file for statistical analysis. This can then be used with anything from Excel to Perseus or a programming environment like R.

Go and see the presentation if you would like to see examples of what you can do with the output. If you would like to know more about crosslinking or the changes to modification detection check out those respective presentations as well.

Protein FDR

- Automatically calculated when decoy option selected
- Displayed in the decoy section:

▼ *Sensitivity and FDR (reversed protein sequences)*

	Target	Decoy	FDR		
Protein family members	4699	219	4.66%		
<input type="text" value="Sequences"/> <input type="text" value="above"/> <input type="text" value="homology"/>	25588	296	1.16%	<input type="button" value="Adjust to"/>	<input type="text" value="1%"/>

Decoy results are available in [the decoy report](#).

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



The new Protein FDR value is automatically calculated when the decoy search option is selected.

It is displayed above the peptide/PSM decoy results.

Based on the following assumptions and definitions:

- Only significant peptide sequence matches (PSMs) are used
- A protein family member may represent multiple same set proteins
- Protein count used for FDR is count of family members
- A protein ID is considered true positive if it contains at least one positive PSM
- A protein ID is a false positive when all the PSM's are false positives

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



Protein FDR is based on the following assumptions and definitions:

By default the protein family report only shows peptide sequence matches (PSMs) with significant scores and are used for the protein family assignment.

A protein family member may represent multiple same set proteins. Only members of all the protein families in the report, those that contain a unique peptide, are counted.

Protein count used for FDR is count of family members. That is, if the report contains 2 families, one with 4 members and the other with a single member, this counts as a total of 5 proteins. Same-set, sub-set and intersection proteins are not counted.

While a protein ID is a false positive when all the PSM's are false positives. Just one true PMS would make the protein identification a true positive. This is very important.

Based on the following assumptions and definitions:

- **Given the number of proteins and the numbers of true and false peptide sequences we use a hypergeometric model to estimate the number of proteins are truly false positive**
 - Simplified approach to that used by MAYU, from the Aebersold group

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



Given the number of proteins and the numbers of true and false peptide sequences we use a hypergeometric model to estimate the number of proteins that are truly false positive

The algorithm is a simplified approach to that used by MAYU, from the Aebersold group

<https://www.mcponline.org/content/8/11/2405><https://www.mcponline.org/content/8/11/2405>

The main differences are that we do not make a separate estimate of the FDR for one-hit wonders and we do not partition the database by protein size. We use a simpler estimate for the number of false proteins in the target database, based on the assumption that the number of decoy proteins never reaches a significant proportion of the database size.

Example

- Target database has 1000 entries
- Search results
 - 500 target proteins
 - 10 decoy proteins
 - FDR $10/500=2\%$?
- False PSMs distributed across 10 proteins
 - Some will also contain true PSMs
 - Half proteins in target database contain true PSMs
 - Estimate that only 5 target proteins contain nothing but false PSMs
- Protein FDR is $5/500 = 1\%$

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



Imagine the and the search results show 500 target proteins and 10 decoy proteins. Does this mean protein FDR is $10/500 = 2\%$? No, it does not. We can assume the false PSMs in the target are distributed across 10 proteins, but some of these will also contain true PSMs, so should not be counted as false. Since half the proteins in the target database contain true PSMs, a reasonable estimate would be that only 5 target proteins containing nothing but false PSMs, so that the protein FDR is $5/500 = 1\%$.

Adjusting the Protein FDR

Default protein FDR

Format	Significance threshold p<	0.05	Max. number of families	AUTO	[help]
	Display non-sig. matches	<input type="checkbox"/>	Min. number of sig. unique sequences	1	
	Preferred taxonomy	All entries	Dendrograms cut at	0	

▼ Sensitivity and FDR (reversed protein sequences)

	Target	Decoy	FDR	
Protein family members	4699	219	4.66%	
Sequences	above	homology	25588	296 1.16% Adjust to 1%

Set Min. number sig. unique peptide sequences to 2

Format	Significance threshold p<	0.05	Max. number of families	AUTO	[help]
	Display non-sig. matches	<input type="checkbox"/>	Min. number of sig. unique sequences	2	
	Preferred taxonomy	All entries	Dendrograms cut at	0	

▼ Sensitivity and FDR (reversed protein sequences)

	Target	Decoy	FDR	
Protein family members	3366	6	0.18%	
Sequences	above	homology	25588	296 1.16% Adjust to 1%

MASCOT : New features in Mascot 2.7

© 2020 Matrix Science



The default significance threshold for a Mascot search is usually 0.05 and this will often give a peptide FDR in the region of 5%. In this dataset the protein FDR is ~4.5%. If we want to lower the Protein FDR we can try adjusting the peptide FDR to 1%. In this case the default peptide FDR is already close to 1% so we can try adjusting the report other ways.

We can adjust the Minimum number of significant unique sequences. This has quite a strong affect on the Protein FDR. If we change it to 2 and eliminate the “one hit wonders” the protein FDR drops to 0.18%.

Adjusting the Protein FDR

Format	Significance threshold p<	0.005	Max. number of families	AUTO	[help]
	Display non-sig. matches	<input type="checkbox"/>	Min. number of sig. unique sequences	1	
			Dendrograms cut at	0	
	Preferred taxonomy	All entries			

▼ Sensitivity and FDR (reversed protein sequences)

	Target	Decoy	FDR				
Protein family members	4198	38	0.91%				
Sequences	above	homology	19483	49	0.25%	Adjust to	1%

Decoy results are available in [the decoy report](#).

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



Alternatively if we can adjust the Significance threshold for the results. I took a guess and reduced it by a factor of 10 from the default values of 0.05 to 0.005 and clicked the format button. The resulting Protein FDR is approximately 1%.

The current HUPO guidelines Interpretation Guidelines for large-scale results recommend adjusting the settings to lower than 1% protein-level global FDR so after formatting this search result would meet those guidelines.

When interpreting the results a protein FDR of 1% only tells us that 1% of the proteins listed are wholly false. This doesn't mean the other 99% are "correct". In particular, where there are same-set proteins, we cannot say which one is "correct". This is because database redundancy causes protein inference ambiguity and we can account for the PSM evidence using several sets of proteins. It is important to remember that a protein accession number in the summary report does not mean "this is the correct protein", it means "the correct protein is likely to be very similar to one of the set of proteins represented by this family member".

Fragment charge states

- Mascot search supports fragment ion charge states of 1+ and 2+
- Top down analysis generates fragment ions with a charge states >2+
- Current solution is to discharge to MH+ for higher charge states

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



Mascot Server automatically searches for fragment ions as singly or doubly charged ions but not higher charge states.

Middle and top down analysis of large peptides and proteins can generate fragment ions with charge states of 3+ or higher.

In the past if you have fragment ions with higher charge states, you needed to deconvolute them to MH+.

This is supported by Mascot Distiller and is still the preferred method.

Fragment charge states

- MGF format now supports fragment ion charge states in the third column

1052.529	265.0801	7	236.4717	7360.649	0.016802	0.710974	232
1068.542	180.6315	7	87.04799	7472.745	0.01753	0.707436	480
1072.044	126.6786	2	87.0784	2142.073	0.019016	0.728941	360
1078.925	274.3603	5	63.94548	5389.586	0.018413	0.74134	238
1082.002	620.7084	9	1.837065	9728.949	0.064662	0.511568	298
1095.742	182.8241	5	117.3447	5473.675	0.018439	0.691048	446
1104.344	446.9772	5	3610.589	5516.682	0.017342	0.954766	116
1111.563	85.08009	1	25.13187	1110.556	0.024389	0.785191	475
1132.535	539.4132	2	1581.047	2263.056	0.018836	0.978732	22
1141.594	91.67654	1	10.8601	1140.587	0.019478	0.885844	338
1141.96	364.0111	5	228.1712	5704.762	0.018316	0.941458	190
1181.524	286.7844	1	340.8262	1180.517	0.020134	0.975948	58
1264.863	256.5222	4	997.5933	5055.423	0.021817	0.950108	378

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



For software that does not decharge the peaks we added support for specifying the charge state in the third column of the peak list. So you have mass, intensity or area and now charge state. Any data in the later columns will be ignored.

Fragment charge states

- New feature on by default in Mascot Server 2.7
- Controlled by “DechargeFragmentPeaks” configuration option
 - This positive integer is the maximum absolute charge state to be decharged
 - Default value “10”

The new feature is on by default and controlled by the “DechargeFragmentPeaks” configuration option.

DechargeFragmentPeaks is the maximum absolute charge state to be decharged and is set to 10 by default.

Better filtering of modifications with complementary deltas in an error tolerant searches

- Up to Mascot Server 2.6 filtering eliminated modifications that are meaningless given the estimated mass error
 - Q->K, in most cases.
- Extended to filter pairs of modifications that sum to $0 \pm$ precursor tolerance
 - Asn->Gly and Gln->Ala are exact inversions of Carbamidomethyl
 - Deoxy, Ser->Ala, and Tyr->Phe are exact inversions of Oxidation

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



The next feature is “Better filtering of modifications with complementary deltas in an error tolerant searches”.

Error tolerant search results automatically remove matches where the delta of the modifications is less than smaller of the precursor or fragment ion tolerances. For example Q->K, in most cases.

In Mascot Server 2.7 we extend that filtering to remove matches with pairs of modifications that where the deltas sum to zero within the precursor tolerance.

Such matches allow the fragment masses to be shuffled around to get a better score while leaving the parent mass unchanged. These are almost always false, with the true match, without either mod, getting a slightly lower score.

Some examples of this are issues are:

Asn->Gly and Gln->Ala are exact inversions of Carbamidomethyl

Deoxy, Ser->Ala, and Tyr->Phe are exact inversions of Oxidation

Better filtering of modifications with complementary deltas in an error tolerant searches

- Mascot Server 2.6

```

4612: Scan 6812 (rt=45.8184) [C:\Users\johnn\Downloads\trypsin\Yeast_In-qel_digest_3.raw]
19526 1106.0533 2210.0920 2210.0967 -2.12 0 52 ▼ 1 U R.LGRHNIDVLEGNQFVNAAK.I + Carbamidomethyl (N-term); [-57.0215 at 85]
-2.11 0 52 2 R.LGRHNIDVLEGNQFVNAAK.I + Deamidated (NQ); [-0.9840 at C-term]
-2.11 0 51 3 R.LGRHNIDVLEGNQFVNAAK.I + Deamidated (NQ); [-0.9840 at C-term]
-2.11 0 50 4 R.LGRHNIDVLEGNQFVNAAK.I + Deamidated (NQ); [-0.9840 at C-term]
-2.12 0 46 5 R.LGRHNIDVLEGNQFVNAAK.I + Carbamidomethyl (N-term); Deamidated (NQ); [-58.0055 at 83]
-2.12 0 45 0.00064 6 R.LGRHNIDVLEGNQFVNAAK.I
4.32 0 2 12 7 IWPNDLLMTEYLCKK + Carbamidomethyl (N-term); Deamidated (NQ)
-0.086 2 0 18 8 LMSITTEQQGQLVAMQREK + 2 Oxidation (M); 3 Deamidated (NQ)
-0.086 2 0 18 8 LMSITTEQQGQLVAMQREK + 2 Oxidation (M); 3 Deamidated (NQ)
-0.086 2 0 18 8 LMSITTEQQGQLVAMQREK + 2 Oxidation (M); 3 Deamidated (NQ)
  
```

- Mascot Server 2.7

```

3597: Scan 5738 (rt=41.3735) [C:\Users\johnn\Downloads\trypsin\Yeast_In-qel_digest_2.raw]
Score > 25 indicates identity
Score > 15 indicates homology
19526 737.7046 2210.0919 2210.0967 -2.19 0 (45) 6.3e-05 ▼ 1 U R.LGRHNIDVLEGNQFVNAAK.I
No other peptide matches in query
  
```

MASCOT : New features in Mascot 2.7

© 2020 Matrix Science



Here I show the reporting of a peptide With a Carbamidomethyl modification at the N-terminal and an equivalent loss of 57Da Asn->Gly modification at the Ans 5 position. Expanding the top 10 matches of the peptide shows the unmodified peptide with a significant score 6th on the list.

In Mascot Server 2.7 all of the top 5 matches have been filtered out, and the bottom 4 too, leaving just the unmodified sequence match.

Combine searches in Mascot Daemon

- IIS web server has a 4GB peak list upload limit.
- “Merge MS/MS files into single search” option can easily exceed this limit.
- There are a number of ways to combine searches after they have been run
 - See How do I do that? - Tips and tricks from an expert from the ASMS 2015 user meeting.
 - http://www.matrixscience.com/workshop_2015.html

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science




We added a new feature to Mascot Daemon which allows you to easily combine search results.

When working with fractionated samples you will want to combine the searches from multiple data files into a single result. The easy way to do this is with the “Merge MS/MS files into single search” option in Mascot Daemon. However if you are using the Windows operating system for your Mascot Server and the default IIS web server, you can hit a 4 GB upload limit.

There are a number of ways around this limit such as using Mascot Daemon on the same computer as Mascot Server, so the searches are submitted on the command line, moving from the IIS web server to Apache or using Mascot Distiller for peak picking and merging the results there. These are described in an older talk from 2015.

Combine searches in Mascot Daemon



MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



Now Mascot Daemon 2.7 gives you another way to merge searches. This is useful for when users are using msconvert or other peak picking third party tool and not using the merging features of Mascot Distiller.

Instead you can select multiple searches in a Mascot Daemon task by CTRL+click individually searches or shift+click a range then right click and choose combined report.

The combined search will open in a web page and list the results files that have been merged at the top of the report.

This will work with searches that have been processed by any peak picking software including Mascot Distiller.

MGF exporting spectra in the original order

The screenshot displays the 'MASCOT Search Results' interface. At the top, there is a header with the Matrix Science logo and the title 'MASCOT Search Results'. Below this, search details are listed, including the user, search title, MS data file, databases, taxonomy, and timestamp. A search bar is visible with filters for 'All', 'Non-significant', and 'Unassigned'. The 'Export' button is highlighted with a red circle. Below the search bar, there are expandable sections for 'Search parameters', 'Score distribution', 'Modification statistics for all protein families', and 'Legend'. The 'Export search results' section is expanded, showing 'Export format' set to 'MGF Peak List' and 'Output spectra in' set to 'original order (for use with mzIdentML or mzTab)'. The 'original order' option is selected with a radio button and highlighted with a red circle. Other options include 'Mascot query order'. At the bottom of this section, there are buttons for 'Show command line arguments' and 'Export search results'. The footer of the page contains the text 'MASCOT : New features in Mascot 2.7', the copyright notice '© 2020 Matrix Science', and the Matrix Science logo.

This next feature requires some background. When a peak list is submitted to Mascot, the neutral molecular mass (M_r) is calculated for each precursor and queries are sorted lowest to highest. The precursors are given a new index number based on this sorting.

You can run into a problem related to the change in order of the precursors in an MGF file when using search results with third party analysis software. When generating a mzIdentML file from the Mascot Server search results the spectrum IDs field in the mzIdentML file refer to the original MGF file precursor order.

Some third-party analysis software, such as the PRIDE submission tool, maps the data in the MGF file to the data in the mzIdentML file. If a collaborator exports the search results from Mascot but does not have access to the original mgf file they can export it from the search results.

In Mascot Server 2.6 the newly created MGF file is exported in the sorted order of the precursors, not the original order, the mapping between the mzIdentML file and the MGF file will fail and the collaborator is left thinking that the either Mascot Server or the third party tool is broken. To prevent this misunderstanding the MGF file can now be exported with the precursors in the original order to maintain mapping compatibility.

To put the new feature to use go to the Export menu, choose MGF Peak list, and it is the default option.

Handle non-standard modification names in a spectral library

- Proteomics Standards Initiative modification names guidelines introduced in 2008.
- Old names or acronyms did not meet the new guidelines.
- Older data sets mined for spectral libraries contained these obsolete names and the spectra were not recognized by Mascot.

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



In 2008 the Proteomics Standards Initiative released community standards for representation of protein modification data by use of a consensus nomenclature. Prior to that there were some differences in modification names between search engines and research groups. Also some of the old names or acronyms did not meet the new guidelines. Some of the modifications had a common name but it did not conform to the new standards, for example "Sodiated" = "Cation:Na⁺". Or used a lab acronym "CAM" = "Carbamidomethyl".

These names live on in old data sets some of which have been used to create spectral libraries that are publicly available from NIST. Mascot Server could not match these modification names to the newer ones and would report an error for those spectra, effectively excluding them from the library.

Handle non-standard modification names in a spectral library

- **mascot/config/library_mod_aliases**
 - Preconfigured with 46 common aliases and acronyms

```
# The default list provided is supplied by NIST
#
"AB_old_ICATd0" = "ICAT-D"
"AB_old_ICATd8" = "ICAT-D:2H(8)"
"Acetylation" = "Acetyl"
"CAM" = "Carbamidomethyl"
"Deamidation" = "Deamidated"
"Dehydro" = "Dehydrated"
"di-Methylation" = "Dimethyl"
"dihydroxy" = "Dioxidation"
"formylkynurenin" = "Dioxidation"
"HMVK86" = "HMVK"
"hydroxykynurenin" = "Trp->Hydroxykynurenin"
```

MASCOT : *New features in Mascot 2.7*

© 2020 Matrix Science



To fix this we have introduced an aliases file that maps the old names on to the new names.

The aliases file resides in the Mascot/config directory and is preconfigured with 46 common aliases and acronyms.

It is extendible as well.