

Maximising Sensitivity with Percolator

MASCOT

MATRIX
SCIENCE

Terminology

| | | The MS/MS spectrum comes from a peptide sequence in the database | |
|--|-------|--|----------------|
| | | True | False |
| Search reports a match to the correct sequence | True | True positive | False positive |
| | False | False negative | True negative |

False Discovery Rate
= $FP / (FP + TP)$


True Positive Rate
= $TP / (TP + FN)$

False Positive Rate
= $FP / (FP + TN)$

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science




Database searching is a statistical process. Most MS/MS spectra do not encode the complete peptide sequence; there are gaps and ambiguities. Hopefully, most of the time, we are able to report the correct match, a true positive, but not always. If the sequence of the peptide is not in the database, and we report a match below our score or significance threshold, that's also OK, and we have a true negative. The other two quadrants represent failure. A false positive is when we report a significant match to the wrong sequence. A false negative is when we fail to report a match even though the correct sequence is in the database. For real-life datasets, where we cannot be certain that all the correct sequences are present in the database, we don't know whether a failure to get a match to a spectrum is a TN or a FN. When we do a decoy search, we make an estimate of TP and FP, and report a false discovery rate, which is defined as the count of significant matches in the decoy sequences divided by the total count of significant matches in both target and decoy.



The characteristic attributes of any scoring algorithm are sensitivity and specificity. That is, you want as many correct matches as possible, that and as few incorrect matches as possible. This curve, that illustrates the relationship between sensitivity and specificity, is called a ROC curve, which stands for Receiver Operating Characteristic.

This plots true positive rate and false positive rate as a function of a discriminator, such as a score threshold. A good scoring scheme will try to follow the axes, as illustrated by the red curve, pushing its way up into the top left corner. A useless scoring algorithm, that cannot distinguish correct and incorrect matches, would follow the yellow dashed diagonal line.

The origin of the ROC curve has unit specificity, i.e. zero false positives, but also zero true positives. Not a useful place to be. The top right of the ROC curve has unit sensitivity, i.e. 100% true positives, but also 100% false positives, which is equally useless. By setting a significance threshold or a score threshold, you effectively choose where you want to be on the curve.



This is another way to look at it. Even the best scoring scheme cannot fully separate the correct and incorrect matches, as shown here in a schematic way. The score distribution for the correct matches, in green, overlaps that of the incorrect matches, in red. The observed score distribution is the sum of these two curves, in black


When we set a score threshold, we are trying to separate the green and red curves as cleanly as possible. But, the lower the threshold, the more incorrect matches are reported. The higher the threshold, the fewer correct matches

But, what if we could find ways to pull these two distributions further apart, or make the distributions narrower? In other words, better resolve the two distributions. This would allow us to improve the sensitivity for a given false discovery rate.

Sensitivity vs. Specificity

Mascot scoring ignores

- Retention time



MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science




This is perfectly possible. There are many observables that the Mascot scoring algorithm doesn't include. For example, HPLC retention time. If the experimental retention times are generally close to the calculated values, we might suspect outliers are false positive matches

Sensitivity vs. Specificity

Mascot scoring ignores

- Systematic mass errors



MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



The more accurate the mass values, the tighter the mass tolerance can be in a Mascot search. But, Mascot only cares about whether the mass values fall within the specified window. In this example, we are searching trap data with a tolerance of ± 0.6 Da. When we look at a strong match, the scatter of fragment mass values appears to be much tighter, maybe ± 0.1 Da, assuming the single high value is random match. When we look at a low scoring, random match, the errors are uniformly scattered across the tolerance window. So, if we had a match that was close to threshold, the scatter on the fragment mass values would be an indication as to whether it was a correct match or not.

Sensitivity vs. Specificity

Mascot scoring ignores

- Counts of modifications

1. [CH60_CR10P](#) Mass: 76623 Score: 404 Matches: 8(5) Sequences: 7(5)
60 kDa heat shock protein, mitochondrial OS=Cricetulus griseus GN=HSPD1 PE=2 SV=1
 Check to include this hit in error tolerant search or archive report

| Query | Observed | Mr(expt) | Mr(calcd) | Delta | Miss | Score | Expect | Rank | Unique | Peptide |
|--|-----------|-----------|-----------|---------|------|-------|----------|------|--------|---|
| <input checked="" type="checkbox"/> 11 | 417.1822 | 832.3498 | 832.3828 | -0.0329 | 0 | 45 | 0.04 | 1 | U | K.APGFGDNR.K |
| <input checked="" type="checkbox"/> 15 | 451.2499 | 900.4853 | 900.5280 | -0.0428 | 0 | 52 | 0.018 | 1 | U | K.LSDGVAVLK.V |
| <input checked="" type="checkbox"/> 16 | 456.7806 | 911.5467 | 911.5804 | -0.0337 | 0 | 59 | 0.0011 | 1 | U | K.VGLQVAVK.A |
| <input checked="" type="checkbox"/> 21 | 480.7447 | 959.4748 | 959.5036 | -0.0288 | 0 | 45 | 0.06 | 1 | U | R.VTDALNATR.A |
| <input checked="" type="checkbox"/> 24 | 595.7855 | 1189.5565 | 1189.6013 | -0.0448 | 0 | 57 | 0.0036 | 1 | U | K.DIGNIISDAMK.K + Methyl (DE) |
| <input checked="" type="checkbox"/> 51 | 1057.0537 | 2112.0929 | 2112.1323 | -0.0394 | 0 | 116 | 4.4e-009 | 1 | U | R.ALMQGVLLDAVAVTNGPK.G |
| <input checked="" type="checkbox"/> 61 | 828.1238 | 2481.3495 | 2481.3942 | -0.0447 | 0 | (12) | 1e+002 | 1 | U | R.TALLDAAGVASLLTTAEAVVTEIPK.E + 2 Methyl (DE) |
| <input checked="" type="checkbox"/> 62 | 828.1322 | 2481.3748 | 2481.3942 | -0.0194 | 0 | 31 | 1.2 | 1 | U | R.TALLDAAGVASLLTTAEAVVTEIPK.E + 2 Methyl (DE) |

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



Here are some results from a search with 3 variable modifications. If we look at the confident matches, most peptides are unmodified. One carries a single modification and a long peptide carries the same modification at two locations.

Sensitivity vs. Specificity

Mascot scoring ignores

- Counts of modifications

| Query | Observed | Mr(expt) | Mr(calc) | Delta | Miss | Score | Expect | Rank | Unique | Peptide |
|--|-----------|-----------|-----------|---------|------|-------|----------|------|--------|--|
| <input checked="" type="checkbox"/> 52 | 1065.0399 | 2128.0653 | 2128.1365 | -0.0712 | 1 | 3 | 1.6e+003 | 1 | | TSVVGQIHTDSMKLVLR |
| <input checked="" type="checkbox"/> 45 | 640.3281 | 1917.9625 | 1917.9624 | 0.0001 | 1 | 12 | 2.3e+002 | 1 | | EQHDRVLELVSQLLR |
| <input checked="" type="checkbox"/> 49 | 1020.9879 | 2039.9613 | 2039.9220 | 0.0392 | 0 | 4 | 1.4e+003 | 1 | | EEEDLLNELDEMEK + 8 Methyl (DE); Carbamidomethyl (N-term) |
| <input checked="" type="checkbox"/> 39 | 752.8643 | 1503.7141 | 1503.7715 | -0.0574 | 1 | 18 | | 79 | 1 | ESRDVLLCHTK + 3 Methyl (DE); Carbamidomethyl (N-term) |
| <input checked="" type="checkbox"/> 34 | 714.8884 | 1427.7623 | 1427.7619 | 0.0004 | 0 | 20 | | 46 | 1 | EEQEAALSAAK + 3 Methyl (DE); Carbamidomethyl (N-term) |
| <input checked="" type="checkbox"/> 40 | 760.8461 | 1519.6777 | 1519.7646 | -0.0869 | 1 | 18 | | 69 | 1 | EKYTLDFTER + 2 Methyl (DE); Carbamidomethyl (N-term) |
| <input checked="" type="checkbox"/> 60 | 789.1094 | 2364.3063 | 2364.2916 | 0.0147 | 1 | 31 | | 1.0 | 1 | LGMIVDALESIRETLEVK + 3 Methyl (DE); Carbamidomethyl (N-term) |

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



Now look down at the low scoring, random matches on the unassigned list. Some are unmodified, of course, but others are heavily modified. One has 8 methyls plus another modification at the terminus. This is to be expected. Peptides that have a large number of potential modification sites support many possible arrangements and permutations of modifications, some of which match quite well by chance. In other words, there are more degrees of freedom. So, if two matches had the same score, and both had 8 Ds and Es, but one was unmodified and the other had 4 methylations, we might feel greater confidence in the match to the unmodified peptide.

Sensitivity vs. Specificity



Anal. Chem. 2002, 74, 5383–5392

Empirical Statistical Model To Estimate the Accuracy of Peptide Identifications Made by MS/MS and Database Search

Andrew Keller,^{*†} Alexey I. Nesvizhskii,^{*†} Eugene Kolker, and Ruedi Aebersold

Institute for Systems Biology, 1441 North 34th Street, Seattle, Washington 98103

Peptide Prophet

- Expectation maximization
- No-enzyme search
- Positive training set: fully tryptic matches
- Negative training set: non-specific matches

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science




The common factor in these properties is that you have to learn how to use them by looking at a set of results of reasonable size, because the rules are likely to change from search to search. Using a count of modifications might not be such a good idea if you were analysing highly modified histones.

The pioneer of using machine learning on a collection of characteristics was Peptide Prophet from the Institute for Systems Biology. This was, and still is, popular for transforming Sequest scores into probabilities.

It takes information about the matches in addition to the score, and uses an algorithm called expectation maximization to learn what distinguishes correct from incorrect matches. Originally, a widely used approach was to run the Sequest search without enzyme specificity and then assume that matches to fully tryptic peptides were correct and matches to non-specific peptides were incorrect.

Sensitivity vs. Specificity



Semi-supervised learning for peptide identification from shotgun proteomics datasets

Lukas Käll¹, Jesse D Canterbury¹, Jason Weston², William Stafford Noble^{1,3} & Michael J MacCoss¹

NATURE METHODS | VOL.4 NO.11 | NOVEMBER 2007 | 923

Percolator


- Support vector machine
- Target decoy search
- Positive training set: high scoring matches from target
- Negative training set: matches from decoy

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



A more recent development has been to use the matches from a decoy database as negative examples for the classifier. Percolator trains a machine learning algorithm called a support vector machine to discriminate between a sub-set of the high-scoring matches from the target database, assumed correct, and the matches from the decoy database, assumed incorrect. Percolator was developed by the MacCoss group at U. Washington. Lukas Kall is now in Sweden, at the University of Stockholm.

Sensitivity optimisation



M. Brosch, L. Yu, T. Hubbard, J. Choudhary, *J Proteome Res* (2009).

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science

**MATRIX
SCIENCE**

This can give very substantial improvements in sensitivity. The original Percolator was implemented mainly with Sequest in mind, but Markus Brosch at the Sanger Centre wrote a wrapper that allowed it to be used with Mascot results and published results such as this. The black trace is the sensitivity using the Mascot homology threshold (MHT) and the red trace is the sensitivity after processing through Percolator (MP). It doesn't work for every single data set. But, when it does work, the improvements can be most impressive.

Those of you who attended this meeting last year will remember that Markus gave a presentation on this topic

(PSM = peptide sequence match, MIT = Mascot identity threshold)

Percolator

The screenshot shows the Mascot MS/MS Ions Search web interface. The 'Decoy' checkbox is highlighted with a red circle. The interface includes fields for 'Your name', 'Email', 'Search Title', 'Database' (SwissProt), 'Enzyme' (Trypsin/P), 'Fixed modifications' (3-Nitrotyrosine (Y), Acetyl (K), Acetyl (N-term), Acetyl (Protein N-term), Amidated (C-term)), 'Variable modifications' (3-Nitrotyrosine (Y), Acetyl (K), Acetyl (N-term), Acetyl (Protein N-term), Amidated (C-term)), 'Peptide tol.' (0.1 Da), 'MS/MS tol.' (0.5 Da), 'Peptide charge' (1+), 'Data format' (Mascot generic), 'Instrument' (ESI-TRAP), and 'Decoy' (checked). The 'Start Search ...' button is visible at the bottom.

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



Using a decoy database is particularly convenient with Mascot, because it can be done automatically as part of any search

Sensitivity optimisation

The screenshot shows the Mascot search results interface. In the 'Protein Family Summary' section, the 'Show Percolator scores' checkbox is highlighted with a red circle. The interface also shows a 'Filter' section with 'Significance threshold p.c. 0.05' and 'Max. number of families: AUTO'. Below the filter, there are tabs for 'Proteins (482)', 'Quantitation (513)', and 'Unassigned (30250)'. The main content area displays 'Protein families 1-10 (out of 482)' with a tree view and a table of protein families.

| Accession | Protein Name | Score | Description |
|------------------|--------------|-------|---|
| 1::TRY1_BOVIN | TRY1_BOVIN | 1597 | |
| 1 2::CP2CT_MOUSE | | 1307 | Cytochrome P450 2C29 OS=Mus musculus GN=Cyp2c29 PE=... |
| 4 2::CP239_MOUSE | | 293 | Cytochrome P450 2C39 OS=Mus musculus GN=Cyp2c39 PE=... |
| 5 2::CP238_MOUSE | | 202 | Cytochrome P450 2C38 OS=Mus musculus GN=Cyp2c38 PE=... |
| 7 2::CP254_MOUSE | | 535 | Cytochrome P450 2C54 OS=Mus musculus GN=Cyp2c54 PE=... |
| 6 2::CP270_MOUSE | | 69 | Cytochrome P450 2C70 OS=Mus musculus GN=Cyp2c70 PE=... |
| 3 2::CY250_MOUSE | | 382 | Cytochrome P450 2C50 OS=Mus musculus GN=Cyp2c50 PE=... |
| 1 2::GRP78_MOUSE | | 1292 | 78 kDa glucose-regulated protein OS=Mus musculus GN=Hsp... |
| 2 2::HSP70_MOUSE | | 353 | Heat shock cognate 71 kDa protein OS=Mus musculus GN=Hsp... |
| 3 2::HS71L_MOUSE | | 165 | Heat shock 70 kDa protein 1-like OS=Mus musculus GN=Hsp... |

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



The developers of Percolator have kindly agreed to allow us to distribute and install Percolator as part of Mascot 2.3. This option is available for any search that has at least 100 MS/MS spectra and auto-decoy results, but it works best if there are several thousand spectra. To switch to Percolator scores, just check the box and then choose Filter. This is the example search that is linked from the MS/MS Summary report help page

Sensitivity optimisation

▼Decoy search summary

| Peptide matches | in target | in Decoy | FDR |
|--|-----------|----------|-------|
| - above identity threshold | 1484 | 15 | 1.01% |
| - above identity or homology threshold | 1837 | 20 | 1.09% |

Decoy results are available in [the decoy report](#).

▼Decoy search summary

| Peptide matches | in target | in Decoy | FDR |
|--|-----------|----------|-------|
| - above identity threshold | 1985 | 19 | 0.96% |
| - above identity or homology threshold | 1985 | 19 | 0.96% |

Decoy results are available in [the decoy report](#).

| Delta M | Score | Expect | Rank | U | 1 | 2 | 3 | Peptide |
|---------|-------|--------|---------|----|---|---|---|---------------|
| 0.1363 | 0 | 33 | 0.00049 | ▶1 | U | ■ | ■ | R.LIGDAAK.N |
| 0.3020 | 0 | 14 | 0.039 | ▶1 | U | ■ | ■ | K.VQVEYK.G |
| 0.2841 | 0 | 17 | 0.018 | ▶1 | U | ■ | ■ | R.QLDSTLQK.S |
| 0.4581 | 0 | 18 | 0.015 | ▶1 | U | ■ | ■ | K.VLEDSLK.K |
| 0.0517 | 0 | 21 | 0.0087 | ▶1 | U | ■ | ■ | K.VLEDSLK.K |
| 0.2227 | 0 | 25 | 0.0031 | ▶1 | U | ■ | ■ | K.IITINDQMR.L |


Score > 13 indicates identity

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



Using the Mascot homology threshold for a 1% false discovery rate, there are 1837 peptide matches. Re-scoring with Percolator gives a useful increase to 1985 matches. Note that, in general, the scores are lower after switching to Percolator. The value in the expect column is the Posterior Error Probability (PEP) output by Percolator. A Mascot score is calculated from this and there is a single score threshold, which we will continue to call the identity threshold, with a fixed value of 13 ($-10 \log 0.05$). By keeping the score, threshold, and expect value consistent, we aim not to break any third party software that expects to find these values.

Terminology: FDR & PEP



I've stolen this slide from the talk Markus gave last year because it makes the difference between FDR and PEP very clear. The vertical dashed line is our significance threshold, chosen to give an acceptable false discovery rate (FDR or q value). This is the ratio of the areas under the black and red curves, B/A . That is, it is a property of the set of matches, not of an individual match.

For any particular match, the chance of it being incorrect, given its score, is the Posterior Error Probability (PEP). This corresponds to the ratio of the heights b/a , although we cannot measure a and b directly.

The Mechanics

- All binaries installed as part of Mascot 2.3
- Currently shipping Percolator 1.14
- After any suitable search:
 1. ms-createpip.exe runs, reading the result file and creating a Percolator input file (*.pip) containing a list of features for every query
 2. Percolator runs, taking input from the *.pip file and writing output to two output files (*.target.pop, *.decoy.pop)
 3. When a report is generated, Mascot Parser transparently opens the *.pop file as required
 4. If you view a report from an old result file that is suitable for Percolator, the report script automatically triggers the creation of *.pip and *.pop files

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



The architecture of the integration between Mascot and Percolator.

Features are the observables, e.g. retention time, mass error, count of modifications or missed cleavages, etc.

The Mechanics

```
Command Prompt
C:\inetpub\mascot\bin>ms-createtip.exe --help
Version 2.3.0
Usage: ms-createtip.exe [OPTION] -i filename

Options:
-h, --help          display this help page and exit
-f, --features      display list of features defined in mascot.dat and exit
--sessionID <id>  not normally used because this is run from command line
-o output_file     default is filename.pip
-q <#queries>      override minimum number of queries set in mascot.dat
-s <#sequences>    override minimum number of sequences set in mascot.dat
-a <feature>       add a feature to the list specified in mascot.dat
-r <feature>       remove a feature to the list specified in mascot.dat
-c                use cached results
-p <interval>     progress reports of Process:XX every <interval> seconds
--nocache         do not use cached results
--version         display version number and exit

Features
retentionTime    Retention time in seconds if available
dm              Calculated minus observed peptide mass in Da
mScore          Mascot score (always on)
lgDScore        Mascot score minus Mascot score of next best non-isobaric
                peptide hit
mCalc           Calculated Mr
charge          Charge
dMppm           Calculated minus observed peptide mass in ppm
absDM           Absolute value of calculated minus observed peptide mass in
                Da
absDMppm        Absolute value of calculated minus observed peptide mass in
                ppm
isoDM           Calculated minus observed peptide mass, after eliminating
                possible isotope errors up to 2 Da, in Da
isoDMppm        Calculated minus observed peptide mass, after eliminating
                possible isotope errors up to 2 Da, in ppm
mc              Number of missed cleavages (always 0 if no enzyme)
varnods         Number of modified sites divided by number of modifiable
                sites
```

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



Configuration information is in mascot.dat. This controls which features are used , paths to executables and other files, logging levels, etc. There is some documentation in the Mascot Setup & Installation manual. You can also get help by executing ms-createtip.exe and percolator.exe with the argument --help

The Mechanics

Creating the input file can be time consuming for a large result file, but is a one-time operation

Defaults are set in mascot.dat

- Whether to show Mascot scores or Percolator scores when report first loaded
- Whether to use retention time information if available
- Which features to include

Some miscellaneous points

Limitations

'Protein Features' carry some risk and are currently not implemented (Mascot 2.3.00)

- Feature is essentially a count of the number of sequences assigned to the parent protein, normalised to the length of the protein. 'To those that have, shall be given'
- Concern 1: There is no analogy of this grouping in the decoy database
- Concern 2: FDR is no longer a true peptide FDR and could be misinterpreted

Only the top ranking match is re-scored

- Never get re-ranking of peptide matches. Scores and expect values for other ranks are pro-rated

Unlikely to succeed if results contain very few good matches

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



We decided not to implement 'protein features' because of concerns that the results could be misleading. Essentially, there is only one protein feature: a count of the number of sequences assigned to the parent protein, normalised to the length of the protein. In biblical terms, 'To those that have, shall be given'.

There are some complications to this. For example, many peptides are found in multiple proteins, so which is the 'true' parent? The longest or the shortest or some average. Normalisation is critical if we want to avoid the 'titan effect', where the very largest proteins are promoted because they randomly match a huge number of peptides.

Another concern was that we may get artefacts because the whole concept of target-decoy validation is peptide-centric. Each peptide sequence match being independent of any other. If you increase the score of a weak match simply because it is found in a protein for which there is strong evidence, the FDR cannot be compared with a conventional, 'pure peptide' FDR

Only the top ranking match to each spectrum is used by Percolator. We tried to include all the significant matches, but couldn't get the stats to work properly. This is something Lukas and colleagues are working on, because there would be a real benefit from allowing Percolator to re-rank matches. For example, the features associated with the rank 1 match might indicate that it is unsafe and should be given a high PEP while the rank 2 match looks great and would get a very low PEP. At present, this change in order cannot happen. If the rank 1 match is given a high PEP then the rank 2 match can only be higher

Finally, you must have a population of good, strong matches to provide a positive training set for the SVM. The larger the data set, the more matches you need.

Limitations



Filter Significance threshold p< Max. number of families [\[help\]](#)
Ions score or expect cut-off
Show Percolator scores [\[help\]](#)

▼Decoy search summary

| Peptide matches | in SwissProt | in Decoy | FDR |
|--|--------------|----------|-------|
| - above identity threshold | 13 | 0 | 0.00% |
| - above identity or homology threshold | 45 | 2 | 4.44% |

Decoy results are available in [the decoy report](#).

Filter Significance threshold p< Max. number of families [\[help\]](#)
Ions score or expect cut-off
Show Percolator scores [\[help\]](#)

▼Decoy search summary

| Peptide matches | in SwissProt | in Decoy | FDR |
|--|--------------|----------|-------|
| - above identity threshold | 25 | 2 | 8.00% |
| - above identity or homology threshold | 25 | 2 | 8.00% |

Decoy results are available in [the decoy report](#).

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



So, for example, if we take the famous T. Rex dataset, where there are only a tiny number of high confidence matches in 48,216 spectra, we don't see any sensitivity improvement. There simply aren't enough good matches for the SVM get traction. But, this is the exception. For a more typical search result, Percolator will give sensitivity a significant boost

Retention Time

RT must be included in the MGF peak list

```
scans=44895  
rtinseconds=4696.366
```

Percolator

1. learns how to predict retention time from the sequences in the search result
2. uses the absolute value of the difference between calculated and observed retention time as a predictive feature

Increases processing time

Can be turned on as default in mascot.dat

```
PercolatorUseRT 1
```

Or, can turn on for individual searches with URL argument

```
percolate_rt=1
```

To use retention time as a feature, the experimental RT values must be present in the MGF peak list. Some peak picking utilities simply embed the RT and scan information as free text in the scan title, which won't work. Percolator fits calculated values to the experimental retention times and then uses the deviations for individual matches as a predictive feature. This increases processing time for Percolator, so it is turned off by default. You can enable it as a global default in mascot.dat, or use a URL argument to enable it for an individual search

Retention Time

Original Mascot
results

▼*Decoy search summary*

| Peptide matches | in SwissProt | in Decoy | FDR |
|--|--------------|----------|-------|
| - above identity threshold | 6198 | 61 | 0.98% |
| - above identity or homology threshold | 6198 | 61 | 0.98% |

Decoy results are available in [the decoy report](#).

After Percolator,
no RT

▼*Decoy search summary*

| Peptide matches | in SwissProt | in Decoy | FDR |
|--|--------------|----------|-------|
| - above identity threshold | 7207 | 74 | 1.03% |
| - above identity or homology threshold | 7207 | 74 | 1.03% |

Decoy results are available in [the decoy report](#).

After Percolator,
with RT

▼*Decoy search summary*

| Peptide matches | in SwissProt | in Decoy | FDR |
|--|--------------|----------|-------|
| - above identity threshold | 7473 | 78 | 1.04% |
| - above identity or homology threshold | 7473 | 78 | 1.04% |

Decoy results are available in [the decoy report](#).

MASCOT : Maximising Sensitivity with Percolator © 2010 Matrix Science



Here is an example where enabling retention time as a feature gives a further useful improvement in sensitivity