

Spectral Library Searching

MASCOT

: *Spectral library searching*

© 2017-2023 Matrix Science



Spectral libraries

- Spectral library contains annotated MS/MS spectra
- Match observed spectra directly to library spectra
- Advantages
 - Faster and more specific than database search
 - Easily search non-tryptic peptides or uncommon modifications
- Disadvantages
 - Only identifies peptides that exist in the library
 - Requires good measurement reproducibility
 - Creating high-quality libraries is time consuming

A spectral library is a collection of annotated MS/MS spectra of peptides. Instead of searching observed spectra against a protein sequence database, you search the observed spectra directly against a spectral library. The observed peaks are compared to the annotated library peaks, then scored in some way based on similarity. Typically, the similarity score takes advantage of peak intensity patterns as well as peak masses. It may also utilise peak annotations, such as giving a higher score to b and y ion matches.

There are several advantages compared to protein sequence databases. A library search is often much faster than a database search, as spectral library typically has orders of magnitude fewer peptides than a tryptic digest of a sequence database. A library search can also be more specific than a database search, for example if it contains previously identified non-tryptic peptides or uncommon variable modifications. Selecting a semi-specific enzyme or including many uncommon variable modifications in a database search greatly increases the search duration. Searching a pre-prepared library of spectra of semi-specific peptides is much faster.

There are no free lunches, so of course spectral library searching has weaknesses. You can only identify peptides for which a spectrum exists in the library. If the library contains a peptide sequence with one missed cleavage, you will not get a match to a peptide with two missed cleavages. A library search also requires decent measurement reproducibility. If peak intensities vary wildly between repeat runs, it is harder to get a good library match. Finally, creating high-quality libraries is time consuming and

some care is needed. Fortunately, Mascot ships with several predefined spectral libraries to help you get started.

MASCOT MS/MS Ions Search

Your name: Email:

Search title:

Database(s):

Taxonomy:

Test_DNA
ZEST_human
Spectral library (SL)
MassIVE_HumanHCD
NIST_HSA_IonTrap
NIST_Human_HCD
NIST_Human_IonTrap
NIST_Rat_QTof
PM_PIZ_SL
PRIDE_Contaminants

Peptide tol. \pm ppm $\#^{13}\text{C}$

Peptide charge:

Data file: No file selected.

Data format:

Instrument:

Decoy: ☐

MS/MS tol. \pm Da

Monoisotopic ☒ Average ☐

Precursor: m/z

Error tolerant: ☐

Target PSM FDR:

MASCOT

: Spectral library searching

© 2017-2023 Matrix Science

**MATRIX
SCIENCE**

Mascot Server can search spectral libraries using MSPepSearch from Steven Stein's group at NIST. When submitting a search, any combination of amino acid FASTA or nucleic acid FASTA databases, and spectral libraries can be selected. Here, we perform a simple search of some data from CPTAC study 6 against a NIST yeast library

Most search parameters – modifications, enzyme, missed cleavages, taxonomy, and instrument – simply don't apply to a library search. All that matters is how well the experimental spectrum matches the one in the library. The main exceptions are the precursor and fragment mass tolerances.

MASCOT Search Results

User : Shirley Jackson
 E-mail :
 Search title : Yeast SL example
 MS data file : klc_031308p_cptac_study6_68011.mgf
 Database : NIST_5.cerevisiae_tonTrap 20180806 (92,609 library entries)
 Timestamp : 13 Jul 2023 at 19:10:42 GMT

Re-search ☒ All ☐ Non-significant ☐ Unassigned [\[help\]](#) Export As XML

Search parameters
 Score distribution
 Modification statistics for all protein families
 Legend

Protein Family Summary

Format Significance threshold 300 Max. number of families AUTO [\[help\]](#)
 Display non-sig. matches ☐ Min. number of sig. unique sequences 1
 Dendrograms cut at 0

Sensitivity

Proteins (976) [Report Builder](#) [Unassigned \(2884\)](#) [6 permalinks](#)

Protein families 1-10 (out of 929)

10 per page 1 2 3 4 5 6 ... 93 [Next](#) [Expand all](#) [Collapse all](#)

Accession contains Find Clear

Family	Accession	Description
1	KPYK1_YEAST	24012 Pyruvate kinase 1 OS=Saccharomyces cerevisiae (strain ...)
2	1 G3P3_YEAST	17399 Glyceraldehyde-3-phosphate dehydrogenase 3 OS=Sacc...
	2 G3P2_YEAST	14945 Glyceraldehyde-3-phosphate dehydrogenase 2 OS=Sacc...
	3 G3P1_YEAST	7854 Glyceraldehyde-3-phosphate dehydrogenase 1 OS=Sacc...
3	1 HSP72_YEAST	16177 Heat shock protein 72 OS=Saccharomyces cerevisiae (...)
	4 BID_YEAST	3054 Bcl-2 interacting domain 1 OS=Saccharomyces cerevisiae (...)

MASCOT : Spectral library searching © 2017-2023 Matrix Science **MATRIX SCIENCE**

On completion of the search, the matches are reported in a protein family summary. In order to generate such a report, we need reliable and accurate protein inference.

Protein inference for library matches

- Library entries are peptides, not proteins, which means that protein information is only present as annotations
- Such annotations are optional, and may be missing
- Even when accessions are present
 - Reliability is unknown
 - Accession may not have any external meaning
 - Will rarely extend to more than a single accession per library entry

There are some difficulties associated with Protein inference for library matches. First of all, library entries are peptides, not proteins, which means that protein information is only ever present as annotations. Such annotations are optional, and may be missing, as in the case of most PRIDE libraries.

Even when present, the reliability is unknown. The accession could be a meaningless number or string. And, I've never seen a library with more than a single accession per library entry, so protein inference will be inaccurate for shared peptides.

Protein inference for library matches

- A reference FASTA database must be specified for each library file
 - Library entries mapped to all proteins in reference that contain the sequence
- If library entry not found in reference database, the accession in the library annotations is used
- If no accession, the peptide sequence is used as the accession

Our solution is to require a reference FASTA database to be assigned to each library file when it is added to the system. The default is SwissProt, with an appropriate taxonomy filter, but any online FASTA database can be chosen. This allows Mascot to map most of the library peptides to accessions in the reference database. This mapping is done at the sequence level, with no constraints from enzyme specificity. If a library entry has a novel sequence, not found in the reference database, the accession in the library annotations is used. If there is no accession, the peptide sequence is treated as the accession, so that duplicate matches to the same peptide can be grouped, if nothing else.

MASCOT Search Results

User : Shirley Jackson
E-mail :
Search title : Yeast SL example
MS data file : klc_031308p_cptac_study6_68011.mgf
Database : NIST_5.cerevisiae_tonTrap 20180806 (92,609 library entries)
Timestamp : 13 Jul 2023 at 19:10:42 GMT

☒ All
 ☐ Non-significant
 ☐ Unassigned
 [\[help\]](#)
 As XML

Search parameters
Score distribution
Modification statistics for all protein families
Legend

Protein Family Summary

Significance threshold
Max. number of families

☐ Display non-sig. matches
 ☐ Min. number of sig. unique sequences

Dendrograms cut at

[\[help\]](#)

Sensitivity

Proteins (976)

[\[permalink\]](#)

Protein families 1-10 (out of 929)

10 per page
 1 2 3 4 5 6 ... 93

1 **KPYK1_YEAST** 24012 Pyruvate kinase 1 OS=Saccharomyces cerevisiae (strain ...

2

1 **G3P3_YEAST** 17399 Glyceraldehyde-3-phosphate dehydrogenase 3 OS=Sacc...
2 **G3P2_YEAST** 14945 Glyceraldehyde-3-phosphate dehydrogenase 2 OS=Sacc...
3 **G3P1_YEAST** 7854 Glyceraldehyde-3-phosphate dehydrogenase 1 OS=Sacc...

3

1 **HSP72_YEAST** 16177 Heat shock protein 72 OS=Saccharomyces cerevisiae (...
4 **BID_YEAST** 3054 Bcl-2 interacting domain-containing protein 1 OS=Saccharomyces...

MASCOT

: Spectral library searching

© 2017-2023 Matrix Science

MATRIX SCIENCE

Here's the library search report again. Protein inference allows us to create a report for library matches that is near identical to a report for FASTA database matches.

Yeast SL example (Mascot Search X)

localhost/mascot/cgi/master_results_2.pl?file=-%2Fdata%2F20230713%2F017

Accession contains Find Clear

	Score	Mass	Matches	Sequences	
2.1 dG3P1_YEAST	17399	35724	49 (49)	28 (28)	Glyceridehydride-3-phosphate dehydrogenase 3 OS=Saccharomyces cerevisiae (strai...
2.2 dG3P2_YEAST	14945	35824	42 (42)	25 (25)	Glyceridehydride-3-phosphate dehydrogenase 2 OS=Saccharomyces cerevisiae (strai...
2.3 dG3P1_YEAST	7854	35728	25 (25)	18 (18)	Glyceridehydride-3-phosphate dehydrogenase 1 OS=Saccharomyces cerevisiae (strai...

Redisplay All None

58 peptide matches (52 non-duplicate, 6 duplicate)

Auto-fit to window

Query Dupes	Observed	Mr (expt)	Mr (calc)	ppm	M	Score	Source	Expect	Rank	U	1	2	3	Peptide
d461	356.193320	710.372088	710.371155	1.31	0	881	SL	7.7e-08	1					K.RIDAGAK.K
d440	386.722490	771.430428	771.431152	-0.94	0	812	SL	3.8e-07	1					N.CLAPLAK.V + Carb...
d493	391.231490	780.448428	780.449356	-1.19	0	795	SL	5.6e-07	1	U				K.RIIVDQK.K
d588	398.211920	794.409288	794.410950	-2.09	0	559	SL	0.00013	1					K.LTGMAPR.V
d511	406.209590	810.404608	810.405746	-1.40	0	451	SL	1.5e-05	1					K.LTGMAPR.V + Oxid...
d502	418.198710	834.363968	834.361954	1.10	0	699	SL	6.4e-06	1					K.YDSTRIR.V
d590	420.714160	839.413768	839.413707	0.013	0	609	SL	4.1e-05	1					K.VDQPSHK.D
d585	424.247500	846.480448	846.471146	11.0	0	317	SL	0.034	1	U				R.LATQFQR.I
d1115	432.731200	863.447848	863.450150	-2.67	0	360	SL	0.013	1	U				K.LATQFQR.D
d1250	440.729020	879.443488	879.444946	-1.66	0	762	SL	1.2e-06	1					K.LATQFQR.D
d1487	303.855100	908.543472	908.544220	-0.82	0	501	SL	0.00049	1	U				K.RIIVDQK.I
d1582	459.760750	917.506948	917.508347	-1.52	0	699	SL	5.1e-06	1	U				K.RIIVDQK.I
d1584	306.843510	917.508702	917.508324	0.52	0	740	SL	2e-06	1	U				K.RIIVDQK.I
d2475	504.776920	1007.539288	1007.539948	-0.66	0	722	SL	3e-06	1					K.IATYQER.D
d2476	336.854640	1007.540292	1007.540024	0.27	0	808	SL	4.2e-07	1					K.IATYQER.D
d2485	505.284620	1008.558688	1008.560349	-1.43	0	557	SL	0.00013	1					V.VDLYVERIAK.A
d2911	526.756780	1051.499008	1051.500748	-1.65	0	558	SL	0.00013	1					M.PMGVREK.V
d3040	533.745390	1065.476228	1065.475952	0.26	0	340	SL	0.02	1					K.IYNSASCTTN.C + C...
d3230	544.818290	1087.622028	1087.624954	-2.69	0	590	SL	6.3e-05	1					M.VGVATHQFR.I
d3231	363.548720	1087.624332	1087.625031	-0.64	0	498	SL	0.00052	1					M.VGVATHQFR.I
d3289	548.286890	1094.559228	1094.560745	-1.39	0	479	SL	0.00081	1	U				K.KATYDQIK.A
d3407	554.820710	1107.626868	1107.628754	-1.70	0	885	SL	7.1e-08	1					R.VVDLYVERIAK.A
d3409	370.216590	1107.627942	1107.628723	-0.71	0	629	SL	2.6e-05	1					R.VVDLYVERIAK.A
d3497	560.747290	1119.480028	1119.483353	-2.97	0	910	SL	4e-08	1					R.YAGRYVHDEK.H
d3498	374.167870	1119.481782	1119.483331	-1.29	0	539	SL	0.0002	1					R.YAGRYVHDEK.H
d3564	376.192180	1125.554712	1125.555527	-0.72	0	589	SL	6.4e-05	1					K.EYTVQRIK.V
d3579	387.694620	773.373586	773.374146	-0.72	0	530	SL	0.00025	1	U				K.AAARQPMK.G
d4028	590.338630	1178.663108	1178.665955	-2.42	0	760	SL	1.3e-06	1					R.VVDLYVERIAK.A
d4030	393.895710	1178.665302	1178.665726	-0.36	0	389	SL	0.0064	1					R.VVDLYVERIAK.A

MASCOT

: Spectral library searching

© 2017-2023 Matrix Science



If only libraries are searched, MS PepSearch scores are converted to arbitrary expect values. A score of 300 becomes an expect value of 0.05 and the maximum score of 1000 becomes an expect of 5E-9.

MASCOT MS/MS Ions Search

Your name: ShirleyJackson Email:

Search title: Yeast SL example

Database(s): FungiDB_EST (NA)
NIST_S.cerevisiae_IonTrap (SL)
UP2311_S_cerevisiae (AA)

> <

UniProt_Fig
Uniprot_Pseudomonasaeruginosa
Uniprot_Rice
UniProt_Rust
UP2494_R_norvegicus
UP290289_Malusdomestica
UP29965_C_sabaeus
UP5226_T_rubripes
UP5640_H_sapiens
UP589_M_musculus
UP675_E_coli_K12

Peptide tol. \pm 10 ppm # ^{13}C 0

Peptide charge: 2+

Data file: klc_031308p_cptac_study6_6B011.mgf

Data format: Mascot generic

Instrument: ESI-TRAP

Decoy: ☐

MS/MS tol. \pm 0.6 Da

Monoisotopic: ☒ Average ☐

Error tolerant: ☐

Target PSM FDR: (no target)


Start Search ... Reset Form

MASCOT : Spectral library searching

© 2017-2023 Matrix Science



Let's expand this example and search a mixture of amino acid and nucleic acid databases with a spectral library.



MASCOT : Spectral library searching

© 2017-2023 Matrix Science



Here is the report from the search. We can see all three databases listed at the top of the result report, and each is assigned an index so that we know where each accession comes from. The top hit has an index 3 which corresponds to the UniProt proteome. There are two important differences between this ‘integrated’ report and a library-only report.

Integrated searches (FASTA database + library)

- **Protein inference**

- Library matches are mapped to accessions from the FASTA database
- Reference database accessions or original library annotations only used where this fails

- **Library match scores**

- Take the set of queries where the library and FASTA database matches agree and the Mascot score is significant
- Find scaling factors for library scores in this set such that their mean and standard deviation are the same as Mascot scores
- Assign expect values based on the scaled scores, using the Mascot expect value formula

MASCOT

: *Spectral library searching*

© 2017-2023 Matrix Science



For protein inference, if the peptide sequence can be mapped to one of the FASTA databases being searched, this becomes the preferred accession. The accession from the reference database is only used when this fails.

In an integrated search, we can use the FASTA database matches to create a simple empirical estimate of library score significance. This is achieved by calibrating library scores based on the set of queries where the library and FASTA database searches return the same match and the Mascot score is significant. The shapes of the library and Mascot score distributions in this set are similar and they often have a fairly high correlation. Next, scale these library scores so that they have the same mean and standard deviation as Mascot scores. This produces values on the same scale as Mascot scores. We can now assign expect values to library matches using the same expression as for Mascot matches.

Accession: 3::N1P8V5
 11 same sets of 3::N1P8V5

▼ 54 peptide matches (46 non-duplicate, 8 duplicate)
☒ Auto-fit to window

Query Dups	Observed	Mr(expt)	Mr(calc)	ppm	M	Score	Source	Expect	Rank	U	Peptide
5352	672.8765	1343.7385	1343.7409	-1.73	0	100	AA	1.6e-07	1	U	R.LTSLMVVAGSDLR.R
5426	452.9921	1354.7545	1354.7567	-1.62	0	582	SL	0.00019	1	U	K.TNNPETLVALRK.A
5427	452.9922	1354.7549	1354.7567	-1.36	0	361	SL	0.017	1	U	P.KTNHPETLVALR.K
5702	465.2392	1392.6958	1392.6971	-0.95	0	330	SL	0.033	1	U	K.NGVHIVFASFIR.T + Oxid
6394	500.5497	1498.6272	1498.6298	-1.72	0	462	SL	0.0022	1	U	R.MNFSHGSEYTHK.S
6395	375.6646	1498.6292	1498.6299	-0.48	0	532	SL	0.00052	1	U	R.MNFSHGSEYTHK.S
6405	750.9274	1499.8402	1499.8419	-1.16	0	591	SL	0.00016	1	U	R.LTSLMVVAGSDLR.T
6406	500.9544	1499.8415	1499.8421	-0.42	0	468	SL	0.0019	1	U	R.LTSLMVVAGSDLR.T
6412	751.4172	1500.8199	1500.8235	-2.42	0	641	SL	5.6e-05	1	U	K.YRNPFPILVTR.C + Cach
6415	501.2814	1500.8224	1500.8234	-0.68	0	481	SL	0.0015	1	U	K.YRNPFPILVTR.C + Cach
6501	759.8380	1517.6615	1517.6634	-1.29	0	59	AA	1.7e-05	1	U	K.EPVSMTDDEAR.I
7417	575.6394	1723.8963	1723.8992	-1.71	0	387	SL	0.01	1	U	K.GVNLPPTDVLPALEK.D
7418	862.9557	1723.8969	1723.8992	-1.32	0	24	NA	0.022	1	U	K.GVNLPPTDVLPALEK
7507	874.9940	1747.9734	1747.9720	0.83	0	821	SL	1.4e-06	1	U	R.GDLGIEPAPEVLAVQK.K
7547	880.9423	1759.8701	1759.8741	-2.23	0	48	AA	0.00025	1	U	K.IENQGVNNFDEILK.V
7917	937.0048	1871.9950	1871.9991	-2.21	0	586	SL	0.00017	1	U	K.SEELYPGRPLAIALDTK.G
7918	625.0063	1871.9970	1871.9993	-1.27	0	530	SL	0.00054	1	U	K.SEELYPGRPLAIALDTK.G
8195	667.7041	2000.0905	2000.0942	-1.86	0	463	SL	0.0021	1	U	R.KSEELYPGRPLAIALDTK.G
8196	1001.0531	2000.0916	2000.0942	-1.27	1	46	AA	0.00044	1	U	R.KSEELYPGRPLAIALDTK.G
8197	501.0305	2000.0929	2000.0943	-0.70	0	435	SL	0.0038	1	U	R.KSEELYPGRPLAIALDTK.G
8206	670.0106	2007.0100	2007.0158	-2.89	0	510	SL	0.00082	1	U	K.PTSTTETVAASAAVAEVR
8233	1011.4730	2020.9314	2020.9377	-3.12	0	755	SL	5	1	U	F.VFEKEPVSMTDDEAR.I
8588	621.0749	2480.2705	2480.2759	-2.18	0	549	SL	0.00037	1	U	K.GVNLPPTDVLPALEKDK
8618	870.1036	2607.2890	2607.2849	1.56	0	426	SL	0.0046	1	U	R.NCTPPTSTTETVAASAVAA

Here, the top hit has been expanded. You can see that the top ranking PSMs come from both library and FASTA database. In most cases, the same match is found in two or all three databases, and the listed match is the one with the lowest expect value. An exception can be seen here for query 8233. This peptide is non-specific at the amino terminus and is only found in the library. It will not be matched in the FASTA database because the enzyme for the search was strict trypsin.

Database Manager

Databases (17)

Parse rules (16)

Scheduled updates (0)

Running tasks (1)

Settings

Fasta

Enable predefined definition

Synchronise custom definitions

Create new

Library

Enable predefined definition

Synchronise custom definitions

Create new

Spectral library filters

Databases and spectral libraries

Name	Mode ?	Type ?	Status			Latest task
3UTRrtFull_ncbi	custom	NA	In use		Deactivate	
contaminants	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
cRAP	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
Mus_musculus_GRCm39_genomic	custom	NA	In use		Deactivate	
NCBIprot	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
PRIDE_Contaminants	predefined	SL	In use	Get new files	Deactivate	Update succeeded (view log)
SARS-CoV-2	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
SwissProt	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
TGAtaAGA_ReadthroughProtein2	custom	AA	In use		Deactivate	
UP5640_H_sapiens	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
UP589_M_musculus	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
UP625_E_coli_K12	predefined	AA	In use	Get new files	Deactivate	Downloading (0.0%)
UP6548_A_thaliana	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
UP9136_B_taurus	predefined	AA	In use	Get new files	Deactivate	Update succeeded (view log)
Vertebrates_EST	predefined	NA	In use	Get new files	Deactivate	Update succeeded (view log)
Human_EST	predefined	NA	Offline	Get new files	Activate	Update succeeded (view log)
UniRef100	predefined	AA	Offline	Get new files	Activate	Update succeeded (view log)

Latest predefined definitions files are from Thu Aug 3 06:04:11 2023 (FASTA databases: databases_20230803T100411.xml) and Thu Oct 13 08:27:38 2022 (spectral libraries: libraries_20221013T122738.xml).

Full database status is available on [the database status page](#).

Refresh

MASCOT

: Spectral library searching

© 2017-2023 Matrix Science



Let's turn our attention to administration aspects. Library files in NIST MSP format are handled in Database Manager much the same as Sequence databases in FASTA format. This slide shows the top level screen of Database Manager, with a mixture of FASTA databases and libraries configured for searching. The 'Type' column shows which are AA or NA FASTA and which are spectral library. Most have 'predefined' configuration settings – that is, Matrix Science maintains a master file of configuration settings that is downloaded by Database Manager.

Mascot Database Manager

Database Manager

Databases (134)

Parse rules (48)

Scheduled updates (0)

Running tasks (0)

Settings

Fasta

Enable predefined definition

Synchronise custom definitions

Create new

Library

Enable predefined definition

Synchronise custom definitions

Create new

Spectral library filters

Enable predefined library definition

Predefined library definitions are configuration entries for the most commonly used, publicly available spectral libraries. Configuration and library files for predefined definitions will be automatically kept up to date as long as the Mascot Server machine is connected to the Internet.

Only one instance of each definition can be enabled at any one time, as database and library definition names need to be unique. You can [copy an existing definition](#) to create more than one instance of a predefined definition.


You can also [use a predefined definition as a template](#) when creating a new definition. Such copies will not be kept up to date with the original predefined definition.

Name	Enable	
NIST_BSA_IonTrap	Enable	Could be synchronised with DFCI .
NIST_C.elegans_IonTrap	Enable	Could be synchronised with DFCI .
NIST_Chicken_IonTrap	Enable	Could be synchronised with DFCI .
NIST_D.rerio_IonTrap	Enable	Could be synchronised with DFCI .
NIST_Drosophila_IonTrap	Enable	Could be synchronised with DFCI .
NIST_E.coli_IonTrap	Enable	
NIST_HSA_IonTrap	Enable	Already set up as NIST_HSA_IonTrap .
NIST_Human_HCD	Enable	Already set up as NIST_Human_HCD .
NIST_Human_HCD_2_good	Enable	
NIST_Human_HCD_3_semitryp	Enable	
NIST_Human_HCD_ITRAQ_1	Enable	
NIST_Human_HCD_ITRAQ_2	Enable	
NIST_Human_HCD_ITRAQ_Phospho	Enable	
NIST_Human_IonTrap	Enable	Already set up as NIST_Human_IonTrap .
NIST_Mouse_HCD	Enable	
NIST_Mouse_HCD_ITRAQ	Enable	
NIST_Mouse_HCD_ITRAQ_Phospho	Enable	


MASCOT

: Spectral library searching

© 2017-2023 Matrix Science



To enable a predefined library is a matter of a few mouse clicks.



If the library you want to search is not on the predefined list, you use the 'Create New' Wizard to configure it as a custom database. A particularly interesting case is if you want to create your own library from Mascot search results. This is easily accomplished, as illustrated in the next few slides. Suppose that we are working on aardvark and want to make a custom library for the aardvark proteome. We choose a name and select 'Create from search results'.

The screenshot displays the Mascot Database Manager web application. The browser window shows the URL `t440p-vrk/mascot/x-cgi/db_manager.pl?dbs.new.config-type=new;dbs.new`. The page title is "Create spectral library from search results".

Left Sidebar:

- Database Manager
 - Databases (9)
 - Parse rules (18)
 - Scheduled updates (0)
 - Running tasks (0)
 - Settings
- Fasta
 - Enable predefined definition
 - Synchronise custom definitions
 - Create new
- Library
 - Enable predefined definition
 - Synchronise custom definitions
 - Create new
 - Spectral library filters

Main Content Area:

Library name:
aardvark

Base directory (?)
`/opt/mascot-2.6-dev/sequence`

Library files will be located in the subdirectory `aardvark` of the base directory. The new directory will be created if it does not already exist.

[Previous](#) [Next](#)

Footer:

MASCOT : Spectral library searching © 2017-2023 Matrix Science **MATRIX SCIENCE**

The next screen just gives an opportunity to change the default location for the files.

The screenshot shows the Mascot Database Manager web interface. The browser address bar displays the URL: `t440p-vrk/mascot/x-cg/db_manager.pl?dbs.new_base_path=%2Fopt%2F`. The page title is "Create spectral library from search results".

Database Manager

- Databases (9)
- Parse rules (18)
- Scheduled updates (0)
- Running tasks (0)
- Settings

Fasta

- Enable predefined definition
- Synchronise custom definitions
- Create new

Library

- Enable predefined definition
- Synchronise custom definitions
- Create new
- Spectral library filters

Create spectral library from search results

Library name:
aardvark

Sequence directory
/opt/mascot-2.6-dev/sequence

Reference database
Please choose a reference database. Where possible, protein accessions for peptides in the spectral library will be taken from the specified Fasta file (the reference database). This will make protein inference more reliable and allows a Protein View report to be displayed for a library hit.

NCBIprot

NCBIprot is larger than 5.0 GB. It is not recommended as a reference database.

Taxonomy
If the selected reference database has taxonomy configured, you can optionally choose a taxonomy for reference accessions.

..... Aardvark

MS/MS tolerance
Please enter estimates for the absolute and relative tolerances of the fragment masses in the library. The tolerances in the Mascot search form apply to the data being searched. A library contains experimental spectra, also subject to mass measurement error. It is better to enter values that are too large rather than too small.

0.1 Da
100 ppm

Previous Create

MASCOT

: Spectral library searching

© 2017-2023 Matrix Science

**MATRIX
SCIENCE**

The reference database is used to assign protein accessions to the library entries. Normally, you wouldn't choose NCBIprot because it is such a large and redundant database. But, since SwissProt only contains 10 aardvark entries, we don't have much choice. We must also provide an estimate of suitable MS/MS tolerances for the library contents. If the search results come from multiple instruments, you need to base this on the least accurate of them.

The screenshot displays the Mascot Database Manager web interface. The browser address bar shows the URL: `t440p-vrk/mascot/x-cgi/db_manager.pl?sub=db%3Aaardvark`. The interface is divided into a left sidebar and a main content area.

Left Sidebar:

- Database Manager**
 - Databases (10)
 - Parse rules (18)
 - Scheduled updates (0)
 - Running tasks (0)
 - Settings
- Fasta**
 - Enable predefined definition
 - Synchronise custom definitions
 - Create new
- Library**
 - Enable predefined definition
 - Synchronise custom definitions
 - Create new
 - Spectral library filters

Main Content Area:

Database: aardvark

Copy **Delete**

Name: aardvark

Database type: Spectral library (created from search results)

Database directory: /opt/mascot-2.6-dev/sequence/aardvark/current

Filename pattern: aardvark_*.msp

Create MSP file from search results

Peptide match filters: (none)

Edit filters

Import search results

The spectral library will be created from Mascot search results. Only results files and peptide matches that pass suitable filtering criteria will be included in the library.

Please configure peptide match filters. After that you can add results to the library.

Footer:

MASCOT : Spectral library searching © 2017-2023 Matrix Science **MATRIX SCIENCE**

Peptide match filters are used to select matches for inclusion in the library. We choose ‘Edit filters’.

Peptide match filters for aardvark

The library must have at least one score or expect value filter, typically `expect < 0.01`.

Each individual filter is in a filter group. To add more filters to the group, use the OR button. To add more groups, use the AND button. The peptide match must pass all filter groups to be accepted, but within each group, only one filter needs to succeed.

To remove a filter, leave its value field empty. To remove a filter group, remove all its filters.

Filters are used in two complementary ways:

1. When Database Manager chooses results files to process, only files that might contain suitable peptide matches are included.
2. When Database Manager loops over peptide matches in a results file, only matches that pass the filter are imported to the library.

For example, if you have a filter `DB = SwissProt` and no other DB filters, then only results files that were searched against SwissProt are processed. (Or in a multi-database search, had SwissProt as one of the databases.) When Database Manager loops over its peptide matches, only those that actually come from SwissProt are imported.

Expect value

AND

Score

AND

Taxonomy ☐ is ☐ is not

MASCOT : Spectral library searching © 2017-2023 Matrix Science

There is a lot of flexibility here. This would be a simple filter for PSMs that can be assigned to a specific organism. We only want strong, confident matches in our library, so we require the match to have an expect value less than 0.01 and a score greater than 50. If the set of search results includes duplicate PSMs, only the one with the highest score goes into the library. We choose Save ...

Database Manager

Databases (10)

Parse rules (18)

Scheduled updates (0)

Running tasks (0)

Settings

Fasta

Enable predefined definition

Synchronise custom definitions

Create new

Library

Enable predefined definition

Synchronise custom definitions

Create new

Spectral library filters

Database: aardvark

Copy

Delete

Name

aardvark

Database type

Spectral library (created from search results)

Database directory

/opt/mascot-2.6-dev/sequence/aardvark/current

Filename pattern

aardvark_*.msp

Create MSP file from search results

Peptide match filters

(expect < 0.01 AND score > 50 AND TAXONOMY is included in ".....
Aardvark")

Edit filters

The spectral library will be created from Mascot search results. Only results files and peptide matches that pass suitable filtering criteria will be included in the library.

Import search results

Waiting for t440p-vrk...

MASCOT : Spectral library searching

© 2017-2023 Matrix Science

MATRIX SCIENCE

Which takes us back to the previous page, and we are ready to import search results.

Mascot Database Manager

Import search results in aardvark

Please enter a date range and an optional filepath wildcard. Database Manager will search ("crawl") for results files matching the wildcard and whose last-modified time is within the date range. Peptide matches in these files will be imported in the library if the file has not been processed yet and if the matches pass the filter criteria.

Results file date range
 From midnight (0:00) on to midnight (23:59) on

Filepath wildcard
 The default is to look in the daily directories of the Mascot data directory: ../data/*/*

By default, a results file will be skipped if it has already been imported in the library. You can override this behaviour by ticking the following box. (The other way to force an already imported file to be processed is to change the filter criteria; this will reset the import status of results files.)

☐ Include files already imported

By default, peptide matches are added to the library and existing entries kept. If a new peptide (sequence + mods) with a higher score is found, it will replace the existing entry of the same peptide. If you tick the following box, the entire library contents will be replaced.

☐ Delete existing library contents

MASCOT : Spectral library searching

© 2017-2023 Matrix Science

MATRIX SCIENCE

The only other thing we need to decide is which search result files to crawl. This can be specified as a date range or a wild card file path or some combination of the two. Finally, we add the import task to the queue and the selected files will be crawled as a background task.

You can even schedule automatic updates for such a database, which means that matches can be imported from new result files, created since the last import.


Cleaning out the Litterbox of Proteomic Scientists' Favorite Pet: Optimized Data Analysis Avoiding Trypsin Artifacts

Matthias Schittmayer,^{†,‡,||} Katarina Fritz,^{†,‡,||} Laura Liesinger,^{†,‡} Johannes Griss,[§] and Ruth Birner-Gruenberger^{*,†,‡}

[†]Research Unit Functional Proteomics and Metabolic Pathways, Institute of Pathology, Medical University of Graz, 8010 Graz, Austria


[‡]Omics Center Graz, BioTechMed-Graz, 8010 Graz, Austria

[§]Department of Dermatology, Medical University of Vienna, 1090 Vienna, Austria

 Supporting Information

ABSTRACT: Chemically modified trypsin is a standard reagent in proteomics experiments but is usually not considered in database searches. Modification of trypsin is supposed to protect the protease against autolysis and the resulting loss of activity. Here, we show that modified trypsin is still subject to self-digestion, and, as a result, modified trypsin-derived peptides are present in standard digests. We depict that these peptides commonly lead to false-positive assignments even if native trypsin is considered in the database. Moreover, we present an easily implementable method to include modified trypsin in the database search with a minimal increase in search time and search space while efficiently avoiding these false-positive hits.

KEYWORDS: proteomics, autolysis protected trypsin, database search, search space restriction, misassigned spectra, false positives



MASCOT : Spectral library searching

© 2017-2023 Matrix Science

 **MATRIX
SCIENCE**

Let's look at a practical example of how these new features might be used. This JPR paper reminded us that sequencing grade trypsin is modified by methylation or acetylation of the lysines. Unless these variable modifications are selected in a search, simply including a contaminants database will not be sufficient to catch all trypsin autolysis peptides. The authors suggested a solution based on editing the sequence of trypsin in the FASTA, replacing K with J, and defining J as the mass of dimethylated lysine. This is fine, as far as it goes, but it misses many of the other modifications that are present, not to mention extensive non-specific cleavage.

Creating a trypsin library


- Download data set from PRIDE
- Find “optimal” set of mods with error tolerant searches
- Search with these mods against SwissProt

Type of search	: MS/MS Ion Search
Enzyme	: semiTrypsin
Fixed modifications	: ⚡Carbamidomethyl (C)
Variable modifications	: ⚡Carbamidomethyl (N-term) , ⚡Methyl (K) , ⚡Methyl (N-term) , ⚡Dimethyl (K) , ⚡Dimethyl (N-term) , ⚡Dehydro (C) , ⚡Deamidated (NQ)
Mass values	: Monoisotopic
Protein mass	: Unrestricted
Peptide mass tolerance	: ± 10 ppm
Fragment mass tolerance	: ± 0.5 Da
Max missed cleavages	: 2
Instrument type	: ESI-TRAP
Number of queries	: 26,505

We downloaded the raw files for one of the data sets in this study from PRIDE and tried a variety of error tolerant searches to discover exactly what was present. Based on these results, we chose these search settings. The enzyme specificity was semiTrypsin because peptides show very extensive C-terminal ‘ragged ends’.

Creating a trypsin library

- Large search space, low sensitivity, but many matches to Trypsin
- Import TRYP_PIG matches as new spectral library “Trypsin”



MASCOT : Spectral library searching

© 2017-2023 Matrix Science



This makes the search space very large, but we do get many matches to trypsin and many modified peptides. The search takes a long time and overall sensitivity is not as good as it would be for a simple search with strict trypsin and only one or two variable modifications.

The answer, of course, is to make a library of the trypsin matches and include this in the vanilla search. This is a very powerful option, since it allows any number of modified and non-specific peptides from any number of contaminants to be intercepted with no increase in the search space.

Benchmark small (Mascot Search X)

localhost/mascot/cgi/master_results_2.pl?file=.%2Fdata%2F20221216%2FF001353.dat 110%


Search title : Benchmark small
MS data file : C:\ProgramData\Matrix Science\Mascot Daemon\MGF40 Benchmark small\mascot_daemon_merge.mgf
Databases : 1: SwissProt 2021_04 (565,928 sequences; 204,173,280 residues)
 2: Trypsin 20221216 (113 library entries)
Timestamp : 16 Dec 2022 at 12:43:40 GMT

Re-search ☒ All ☐ Non-significant ☐ Unassigned [\[help\]](#) Export As XML

▼ Search parameters
Type of search : MS/MS Ion Search
Enzyme : Trypsin/P
Fixed modifications : [Carbamidomethyl \(C\)](#)
Variable modifications : [Oxidation \(M\)](#)
Mass values : Monoisotopic
Protein mass : Unrestricted
Peptide mass tolerance : ± 20 ppm
Fragment mass tolerance : ± 0.5 Da
Max missed cleavages : 1
Instrument type : ESI-TRAP
Number of queries : 99,299

► Score distribution
▼ Modification statistics for all protein families

Modification	Delta	Type	Site	Total matches
Carbamidomethyl	57.021464	fixed	C	1357
Oxidation	15.994915	variable	M	505
Deamidated	0.984016	SL	N	127
Dimethyl	28.0313	SL	K	8
Deamidated	0.984016	SL	Q	5
Methyl	14.01565	SL	K	4
Carbamidomethyl	57.021464	SL	E	2
Carbamidomethyl	57.021464	SL	N-term	2
Dehydrated	-18.010565	SL	C	2
Carbamidomethyl	57.021464	SL	C	1

MASCOT : Spectral library searching © 2017-2023 Matrix Science 

Here, we search Swissprot plus the tryptic autolysis library with strict trypsin and a single variable mod - yet still obtain matches to all the modified and non-specific trypsin autolysis peptides.

Benchmark small (Mascot Search: X)

localhost/mascot/cgi/master_results_2.pl?file=.%2Fdata%2F20221216%2F001353.dat% 90%

Accession

contains

Find

Clear

	Score	Mass	Matches	Sequences	emPAI	
1.1	21980	25078	689 (689)	17 (17)	141.88	Trypsin OS=Sus scrofa OX=9823 PE=1 SV=1
1.2	1266	26927	50 (50)	6 (6)	1.90	Trypsin-2 OS=Homo sapiens OX=9606 GN=PRSS2 PE=1 SV=1
1.3	845	27111	37 (37)	5 (5)	1.89	Trypsin-1 OS=Homo sapiens OX=9606 GN=PRSS1 PE=1 SV=1

Redisplay

All

None

751 peptide matches (49 non-duplicate, 702 duplicate)

☒ Auto-fit to window

Query Dups	Observed	Mr (expt)	Mr (calc)	ppm	M	Score	Source	Expect	Rank	U	1	2	3	Peptide
#35618 3	600.3083	1198.6021	1198.6022	-0.15	0	65	AA	7.3e-05	1	U				K.VYNTVDWIK.D
#56785 2	716.8747	1431.7349	1431.7348	0.078	0	241	SL	0.00031	2	U				VLEGNQFVINAAK.I
#64125	763.8844	1525.7542	1525.7435	6.99	0	63	SL	0.027	1	U				K.SGGSTPSPILLQCLK.A + Carbamidomethyl
#64984 2	770.9321	1539.8497	1539.8508	-0.73	0	50	AA	0.00026	1	U				R.VSTISLPTAPPATGK.C
#64992 1	514.2911	1539.8514	1539.8508	0.40	0	26	AA	0.026	1	U				R.VSTISLPTAPPATGK.C
#69486 3	804.4077	1606.8009	1606.8013	-0.30	0	240	SL	0.00031	1	U				N.FNGNTLNDIMLIK.L
#69540 5	804.8999	1607.7853	1607.7854	-0.029	0	318	SL	4.4e-05	1	U				N.FNGNTLNDIMLIK.L + Deamidated
#72330 2	830.3964	1658.7783	1658.7797	-0.83	0	116	AA	4.8e-10	1	U				K.ITSNMFCVFLGGK.D
#72402 3	830.9299	1659.8452	1659.8457	-0.33	0	247	SL	0.00026	1	U				N.IDVLEGNQFVINAAK.I
#73239	838.3940	1674.7735	1674.7746	-0.63	0	16	AA	0.047	1	U				K.ITSNMFCVFLGGK.D + Oxidation (M)
#73805 3	843.9017	1685.7889	1685.7906	-0.98	0	110	AA	2.9e-10	1	U				K.ITSNMFCVFLGGK.D
#74563 3	851.8992	1701.7838	1701.7855	-0.96	0	75	AA	1.3e-06	1	U				K.ITSNMFCVFLGGK.D + Oxidation (M)
#75100 3	857.4069	1712.7993	1712.7995	-0.14	0	245	SL	0.00028	1	U				R.LGHRNIDVLEGNQFVINA.K
#77710 5	887.9516	1773.8886	1773.8886	0.055	0	320	SL	4.2e-05	1	U				H.MIDVLEGNQFVINAAK.I
#78418	897.9101	1793.8056	1793.8065	-0.54	0	47	SL	0.04	1	U				R.SCAAGTECLISGNGNTK.S + 2 Dehydrated; Dimethyl
#83420 11	965.4534	1928.8922	1928.8927	-0.28	0	245	SL	0.00028	1	U				K.IITHNPFNGLNDIM.L + Deamidated
#83739 12	970.9703	1939.9261	1939.9265	-0.22	0	388	SL	7.5e-06	1	U				R.LGHRNIDVLEGNQFVINA.A
#85509 1	1006.4896	2010.9646	2010.9635	0.55	0	214	SL	0.0006	1	U				R.LGHRNIDVLEGNQFVINA.A
#86963 1	1042.0079	2082.0012	2082.0007	0.24	0	304	SL	6.2e-05	1	U				R.LGHRNIDVLEGNQFVINA.K
#87233 1	700.0110	2097.0112	2097.0110	0.081	0	268	SL	0.00015	1	U				G.KHNIDVLEGNQFVINA.K + Carbamidomethyl
#89578 227	737.7044	2210.0915	2210.0967	-2.35	0	100	AA	2.3e-08	1	U				R.LGHRNIDVLEGNQFVINA.K.I
#89832 14	1106.0550	2210.0954	2210.0967	-0.58	0	132	AA	1.6e-11	1	U				R.LGHRNIDVLEGNQFVINA.K.I
#89942 4	553.5313	2210.0960	2210.0967	-0.32	0	40	AA	0.0017	1	U				R.LGHRNIDVLEGNQFVINA.K.I
#90099 4	1106.5484	2211.0822	2211.0795	1.22	0	269	SL	0.00015	1	U				R.LGHRNIDVLEGNQFVINA.K.I + Deamidated
#90102 24	738.0358	2211.0856	2211.0791	2.91	0	319	SL	4.3e-05	1	U				R.LGHRNIDVLEGNQFVINA.K.I + Deamidated

MASCOT

: Spectral library searching

© 2017-2023 Matrix Science

MATRIX SCIENCE

This removes 776 spectra which otherwise might have given rise to false positives.

If you're wondering about the ridiculous emPAI value, it's because the assumption behind emPAI is strict tryptic cleavage. However, the library search is giving all kinds of semitryptic matches, so the model assumptions are not satisfied.

Summary

- Mascot Server uses NIST MSepSearch for spectral library searches
- You can search any combination of FASTA databases and spectral libraries
- Results are presented using the protein family summary report
- A reference FASTA database is assigned to each library file to ensure accurate protein inference
- For an integrated search, library match expect values are determined from the set of matches that have significant Mascot score and where the library and FASTA database searches agree
- MSP files are configured and updated just like FASTA databases
- Libraries can be created by importing results from searches against FASTA databases

MASCOT

: *Spectral library searching*

© 2017-2023 Matrix Science



To summarise.