

Search Parameters

MASCOT

 MATRIX
SCIENCE

Search Parameters

The image displays two side-by-side screenshots of the Mascot search interface. The left screenshot shows the 'MASCOT Peptide Mass Fingerprint' form, and the right screenshot shows the 'MASCOT MS/MS Ions Search' form. Both forms include fields for user name, email, search title, database, enzyme, taxonomy, fixed and variable modifications, protein mass, peptide tolerance, mass values, data file, data format, instrument, and a 'Start Search' button.

MASCOT : Search Parameters

© 2007-2014 Matrix Science



In this presentation, we will describe each of the Mascot search parameters.

If you submit a search from a web browser, you have a choice of three different search forms. All three forms submit to the same search engine, but they have been optimised for three different types of search. The form for a peptide mass fingerprint is shown on the left, and the form for a search of uninterpreted MS/MS data on the right. Most of the controls are common to both.

Search Parameters

The screenshot shows the 'MASCOT Sequence Query' web form. At the top, it has fields for 'Your name' (Liu Scone) and 'Email' (liu@es.edu). Below these are 'Search title' and 'Database(s)' (Invertebrates_EST, Human_EST, Fungi_EST, Environmental_EST, Bacteria). The 'Enzyme' is set to 'Trypsin'. There are checkboxes for 'Allow up to' (missed cleavages) and 'Quantification' (none). Under 'Taxonomy', it says 'All entries'. There are sections for 'Fixed modifications' and 'Variable modifications', both currently set to 'none selected'. A list of modifications is visible on the right, including Acetyl (K), Acetyl (S-term), Acetyl (Protein N-term), Amidated (C-term), Amidated (Protein C-term), Aminoxylation (S-term C), Brn (K), Brn (S-term), Carboxymethyl (C), Carboxymethyl (K), and Carboxymethyl (S-term). The 'Peptide mol. wt.' is 1.2 Da, 'pI' is 4.35, 'MS/MS mol. wt.' is 0.6 Da, and 'Peptide charge' is 10. The 'Instrument' is 'Default', 'Decoy' is checked, and 'Report top' is set to 'AUTO' with 'hits' selected. At the bottom are 'Start Search' and 'Reset Form' buttons.

MASCOT : Search Parameters

© 2007-2014 Matrix Science



The third form is for a sequence query, such as a sequence tag search. The controls on this form are very similar to those on the MS/MS form. The main difference is that we have a text area to type in the queries, rather than a data file upload control.

Help

PMF ✓ SQ ✓ MS/MS ✓

Click on any link for help

MASCOT : Search Parameters

© 2007-2014 Matrix Science

At the top of each slide, there is a key to show which search parameter applies to which type of search.

The labels on the search form are hyperlinks. Just click on them to get detailed help

User details and title

PMF✓ SQ✓ MS/MS✓

Your name	<input type="text" value="Expert User"/>	Email		<input type="text" value="smartie@matrixscience.com"/>
Search title	<input type="text" value="Arabidopsis sample #3476"/>			

- Search form will 'remember' user name and email address in cookie
- If Mascot security is enabled, then this information taken from user database
- Email address used for sending results
- Search title is shown in the report, and can help locate a search in the search log.

MASCOT : Search Parameters

© 2007-2014 Matrix Science

**MATRIX
SCIENCE**

At the top of the form are a couple of fields for user information. The name and email are saved as a browser cookie when a search is submitted, so you don't need to complete them every time.


If you have an in-house server, and Mascot security is enabled, these fields will be populated automatically with the details of the user who is logged in

When you use the Matrix Science public web site, you have to supply a name and email address. This is to allow the results of a search to be returned by email. Usually, search results are returned promptly to your browser window. However, if your connection to the web site is broken before the search is complete, they will be emailed to the supplied address. If you have an in-house server, you can enable this if you wish. It is turned off by default

The search title is free text. You don't have to enter anything. However, it is a good idea to fill in all of these fields, because it makes it much easier to find your old search results in the search log.

Database

PMF✓ SQ✓ MS/MS✓




Choose the right database

- Swiss-Prot good for PMF and MS/MS of well characterised organisms
- NCBIInr or UniRef100 if you want to search all known protein sequences
- ESTs for MS/MS if genome not sequenced

MASCOT : Search Parameters

© 2007-2014 Matrix Science



Choosing the right database is so important that there will be a complete presentation on this topic.

Very briefly, for a peptide mass fingerprint, search a comprehensive, non-redundant database, like SwissProt. If the data are any good, it won't matter if one or two mass values fail to find matches. The advantage of searching a small database is that the search is fast and the reports are concise.

For MS/MS of a well characterised organism, such as human or mouse or yeast, SwissProt is still a good choice. In other cases, search a comprehensive, non-identical database, where every single peptide is explicitly represented, such as NCBIInr or UniRef100.

If the genome of your organism of interest has not been sequenced, it won't be represented in the protein databases, but there may be lots of Expressed Sequence Tags (ESTs). Not advisable for PMF, because many sequences correspond to protein fragments.

In Mascot 2.3 and later, you can select multiple databases for a search. This is particularly useful when you want to search a single organism database and include the sequences of common contaminants, such as BSA and trypsin. One restriction is that you cannot mix AA and DNA databases.

Taxonomy

PMF✓ SQ✓ MS/MS✓

Taxonomy All entries

- Speeds up the search
- Simplifies the result report
- The drop-down list is easily configurable.
- Make sure that the taxonomy indexes are kept up to date.

MASCOT : Search Parameters

© 2007-2014 Matrix Science

MATRIX SCIENCE

If a database contains taxonomy information, we can use this to restrict the search to entries for a particular organism or family. This speeds up the search because, in effect, it makes the database smaller.

Limiting the taxonomy simplifies the result report, because you don't see all the homologous proteins from other species.

The drop down list in the search form is configurable. If you are working on a particular organism, you can easily add this to the list

It is important that the taxonomy is as accurate as possible, which means keeping the indexes up to date

Taxonomy

PMF✓ SQ✓ MS/MS✓

Mascot search status page: https://www.matrixscience.com/x-cgi/ms-status.exe?Show=MAIN_PAGE

Database	File Name	Path Name	Status	State Time	Searches	Map	Unmap	Locked	Threads	Tax IDs	Count
HCBins	HCBins_201407.fasta	/srv/mascot/sequence/HCBins/current/HCBins_201407.fasta	In use	Sun Aug 31 09:53:59	0	YES	NO	YES	-1	4912	15665520
HCBins	HCBins_20140323.fasta	/srv/mascot/sequence/HCBins/current/HCBins_20140323.fasta	In use	Sun Aug 31 09:53:20	1	YES	NO	NO	-1	13091240	6251692
contaminants	contaminants_20090624.fasta	/srv/mascot/sequence/contaminants/current/contaminants_20090624.fasta	In use	Sun Aug 31 09:53:57	0	YES	NO	NO	-1	1011474	593031
cRAP	cRAP_20090731.fasta	/srv/mascot/sequence/cRAP/current/cRAP_20090731.fasta	In use	Sun Aug 31 09:53:57	0	YES	NO	YES	-1	266964	216135
Environmental EST	Environmental_EST_120.fasta	/srv/mascot/sequence/Environmental_EST/current/Environmental_EST_120.fasta	In use	Sun Aug 31 09:53:57	0	YES	NO	YES	-1	131122	112971

MASCOT : Search Parameters

© 2007-2014 Matrix Science



From time to time, its a good idea to go to the database status page and check the stats file for each database. The stats file contains lots of useful information, like whether entries contain illegal characters or whether an entry is too long.

It also tells you how good your taxonomy is. Here are the numbers for the nr database on our web site in March 2014. There are 38 million entries, and 4912 have no taxonomy. In other words, 99.99% of the entries have a taxonomy assigned. If you look at your stats file and see that (say) 10% of the entries have no taxonomy, that's 10% of the entries that are going to be missed whenever you do a search with taxonomy specified.

Taxonomy

PMF✓ SQ✓ MS/MS✓

In most cases, if the correct protein is not in the database, you'd like to see the closest match ... whatever the species

Statistic	Value
Time files compressed	Mon Jul 13 06:51:18 2014
Time files compressed (ms)	140320718
Time / date of fasta file	Mon Jul 9 12:58:00 2014
Time of fasta file (ms)	1404907080
Number of residues	194259949
Number of sequences	546000
Number with invalid residues	0
Number of sequences too long	0
Length of longest sequence	33213
Maximum Accession Length	11
Version of this file	2.4.109
Version of fasta file	6
Type of fasta file	AA
Fasta rule for accession	>.[!"] ' \\ ' ']*
Sequences with invalid taxon names	0
Run sequences for taxonomy	All entries=546000
Run sequences for taxonomy	Archaea (Archaeobacteria)=13270
Run sequences for taxonomy	Eukaryota (Eucaryotes)=17682
Run sequences for taxonomy	Alveolata (Alveolates)=1063
Run sequences for taxonomy	Pinnaeolus (Pinnaeolus)=312
Run sequences for taxonomy	Other Alveolata=751
Run sequences for taxonomy	Metazoa (Metazoa)=103619
Run sequences for taxonomy	Cnidaria (Cnidaria)=3465
Run sequences for taxonomy	Tentaculata (Tentaculata)=5552
Run sequences for taxonomy	Chordata (vertebrates and relatives)=81211
Run sequences for taxonomy	Other vertebrates=83609
Run sequences for taxonomy	Lobe-finned fish and tetrapod clade=76457
Run sequences for taxonomy	Mammalia (mammals)=64263
Run sequences for taxonomy	Primates=42935
Run sequences for taxonomy	Other primates=4725
Run sequences for taxonomy	Rodentia (Rodentia)=24259
Run sequences for taxonomy	Mus=14729
Run sequences for taxonomy	Mus musculus (house mouse)=14678
Run sequences for taxonomy	Rattus=7360
Run sequences for taxonomy	Other rodentia=1610
Run sequences for taxonomy	Other mammalia=13049
Run sequences for taxonomy	Xenopus laevis (African clawed frog)=3385
Run sequences for taxonomy	Other lobe-finned fish and tetrapod clade=8809
Run sequences for taxonomy	Actinopterygii (ray-finned fishes)=3132
Run sequences for taxonomy	Takifugu rubripes (Japanese Pufferfish)=173

MASCOT : Search Parameters

© 2007-2014 Matrix Science



A word of warning. Don't specify a very narrow taxonomy in a search.

Think carefully about what you are trying to achieve when you do this.

If the correct protein from the correct species is not in the database, wouldn't you want to see a good match to a protein from a similar species?

This is especially important for poorly represented species. For example, look at these numbers for the Swiss-Prot 2014_07: half a million entries; 26 thousand entries for rodents, but only 1600 are not either mouse or rat. So, even if you are studying hamster or porcupine, you don't want to choose 'Other rodentia'.

Enzyme

PMF✓
SQ✓
MS/MS✓


Enzyme Trypsin/P

Allow up to 1 missed cleavages

- First choice should normally be the enzyme actually used, and 1 missed cleavage
- Large number of missed cleavages, try increasing to 2
- Use semi-trypsin rather than no enzyme
- No enzyme only in exceptional cases, and never for PMF.

MASCOT : Search Parameters

© 2007-2014 Matrix Science




All the search forms have a drop down list for choosing an enzyme. If your peptides come from an enzyme digest, you need to know what the enzyme was and then choose it from the list.

Setting the number of allowed missed cleavage sites to zero simulates a limit digest. If you are confident that your digest is perfect, with no partial fragments present, this will give maximum discrimination and the highest score for a peptide mass fingerprint.

If experience shows that your digest mixtures usually include some partials, that is, peptides with missed cleavage sites, you should choose a setting of 1, or maybe 2 missed cleavage sites. Don't specify a higher number without good reason, because each additional level of missed cleavages increases the number of calculated peptide masses to be matched against the experimental data. In other words, the missed cleavage parameter should be set by looking at some successful search results to see how complete your digests really are.

Enzyme

PMF ✓ SQ ✓ MS/MS ✓



Olsen, J. V., Ong, S.-E. and Mann, M., *Mol. and Cellular Proteomics*, 3, 608-14 (2004)

MASCOT : Search Parameters

© 2007-2014 Matrix Science



Although some people like to perform searches without enzyme specificity, and then gain confidence that a match is correct if the match is tryptic, this isn't a good idea. If there is evidence for a lot of non-specific cleavage, then a semi-specific enzyme allows one end of the peptide to be non-specific, but not both. Only abandon enzyme specificity if you have no other choice, such as when searching endogenous peptides.

You cannot perform a no-enzyme peptide mass fingerprint. It simply won't work, even if you have good mass accuracy

There is some controversy over the level of non-specific peptides that can be expected in a tryptic digest. Our experience is that the levels of non-specific peptides are very low, less than 3%, unless there is something seriously wrong with the trypsin or the protocol.

Why do we advise so strongly against no-enzyme searches?

Enzyme


PMF✓ SQ✓ MS/MS✓

Fixed modifications: Carbamidomethyl (C)
Variable modifications: Oxidation (M)
Peptide mass tolerance: ± 10 ppm (# 13C = 1)
Fragment mass tolerance: ± 0.1 Da
Max missed cleavages: 2
Instrument type: ESI-TRAP
Number of queries: 44,894
Peptide FDR: 1%

CLE	candidate peptides	seconds	average identity score	matches above identity	Matches above homology
Trypsin	4.4E6	42	26	16,767	17,437
Semi-trypsin	6.9E7	150	38	12,732	15,242
none	3.9E8	670	44	10,681	14,074

MASCOT : Search Parameters

© 2007-2014 Matrix Science

 MATRIX SCIENCE

Here are some numbers for an Orbitrap dataset when we search using strict trypsin, semi-specific trypsin, and no enzyme specificity

As you can see, the no enzyme search takes a lot longer and we get fewer reliable matches.

The reason is simple, the search space for a no-enzyme search is much, much larger than for a tryptic search. This means that the thresholds are higher and we lose marginal matches. Unless you have a high level of non-specific peptides, you lose more than you gain.

So, doing a no-enzyme search in Mascot is not a good idea unless there is a very high level of non-specific peptides. Semi-trypsin will be a better choice if the peptides came from a tryptic digest but there is a high level of non-specific cleavage. Only use no enzyme if the peptides are not the products of an enzyme digest, e.g. MHC peptides or endogenous peptides.

Enzyme

PMF✓ SQ✓ MS/MS✓

Title	Sense	Cleave at	Restrict	Independent	Semispecific	Edit Delete
Trypsin	C-Term	KR	P	no	no	Edit Delete
Trypsin/P	C-Term	KR		no	no	Edit Delete
Arg-C	C-Term	R	P	no	no	Edit Delete
Asp-N	N-Term	BD		no	no	Edit Delete
Asp-N_ambic	N-Term	DE		no	no	Edit Delete
Chymotrypsin	C-Term	FLWY	P	no	no	Edit Delete
CNBr	C-Term	M		no	no	Edit Delete
CNBr+Trypsin	C-Term	KR	P	no	no	Edit Delete
Formic_acid	N-Term	D		no	no	Edit Delete
Lys-C	C-Term	K	P	no	no	Edit Delete
Lys-C/P	C-Term	K		no	no	Edit Delete
LysC+AspN	N-Term	BD	P	no	no	Edit Delete
Lys-N	N-Term	K		no	no	Edit Delete
PepsinA	C-Term	FL		no	no	Edit Delete
semiTrypsin	C-Term	KR	P	no	yes	Edit Delete
TrypChymo	C-Term	FLKRWY	P	no	no	Edit Delete
TrypsinMSIP1	N-Term	J		no	no	Edit Delete
TrypsinMSIP1	C-Term	KR	P	no	no	Edit Delete
TrypsinMSIP1/P	N-Term	J		no	no	Edit Delete
TrypsinMSIP1/P	C-Term	KR		no	no	Edit Delete
V8-DE	C-Term	DEZ	P	no	no	Edit Delete
V8-E	C-Term	EZ	P	no	no	Edit Delete
NoCleave	C-Term	J	ABCDEFGHIKLMNOPQRSTUVWXYZ	no	no	Edit Delete
None						

[Add new enzyme](#) [Main menu](#)

MASCOT : Search Parameters

© 2007-2014 Matrix Science



The list of enzymes is user configurable. Standard entries are described in the help. If you wish, you can modify the definitions or create new ones using the configuration editor.

Mascot supports two categories of mixed enzyme definitions. An independent mixed enzyme is used where multiple sample aliquots have been digested separately, and the digests combined for analysis. This means that the sample could contain (say) tryptic peptides and Asp-N peptides, but no peptides that are tryptic at one end and Asp-N at the other. The second category simulates a single sample aliquot being digested simultaneously or serially by more than one cleavage agent. For example CNBr followed by trypsin.

Remember that enzyme type None simulates cleavage at every peptide bond. For top down searches, where you don't want any cleavage, choose NoCleave.

Enzyme

PMF✓ SQ✓ MS/MS✓

The screenshot shows the MASCOT Sequence Query web interface. The form includes fields for 'Your name' (Liu, Siana), 'Email' (liu@es.ady), and 'Search title' (Glu-fib). The 'Database(s)' dropdown is set to 'Human_EST', 'Pump_EST', 'Environmetal_EST', 'SwissProt', and 'NCBI'. The 'Enzyme' dropdown is set to 'TrypsinP'. The 'Allow up to' dropdown is set to '1 missed cleavages'. The 'Quantification' dropdown is set to 'None'. The 'Taxonomy' dropdown is set to 'Homo sapiens (human)'. The 'Fixed modifications' dropdown is set to 'none selected'. The 'Variable modifications' dropdown is set to 'none selected'. The 'Display all modifications' checkbox is checked. The 'Peptide list, s' dropdown is set to '0.3'. The 'Peptide charge' dropdown is set to '1+'. The 'Peptide mass' dropdown is set to '1369.7'. The 'MS/MS list, s' dropdown is set to '0.5'. The 'Query' text area contains the sequence 'Glu-fib'. The 'Instrument' dropdown is set to 'Default'. The 'Report type' dropdown is set to 'AUTO'. The 'Start Search' button is visible at the bottom.

MASCOT : Search Parameters

© 2007-2014 Matrix Science




Remember that enzyme specificity also applies to Sequence Queries

Quite often, we receive a support email along the lines of "Mascot is broken. I did a search for this peptide and I know its in the database but Mascot failed to find it"

For example, here's a search for glu-fib, a very common sequencing standard. The mass is correct and the sequence is correct. But, when we do a search of SwissProt ...

Enzyme

PMF✓ SQ✓ MS/MS✓



MASCOT : Search Parameters

© 2007-2014 Matrix Science



No results!
Why?

Enzyme

PMF✓ SQ✓ MS/MS✓

The screenshot shows the Mascot Search Results page for the protein FIBB_HUMAN. The page includes the following information:

- Protein View: FIBB_HUMAN**
- Fibrinogen beta chain OS=Homo sapiens GN=FGB PE=1 SV=2**
- Database:** UnimProt
- Score:** 72
- Expect:** 0.0013
- Nominal mass (M₀):** 55902
- Calculated pI:** 5.54
- Taxonomy:** [Homo sapiens](#)
- Sequence similarity is available as an NCBI BLAST search of FIBB_HUMAN against nr.**
- Search parameters**
 - Enzyme:** semiTrypsin: cuts C-term side of KR unless next residue is P. Cleavage is semi-specific. (Peptide can be non-specific at one terminus only.)
- Protein sequence coverage: 2%**
- Matched peptides shown in *bold red*.**
- Unformatted sequence string:** [SLIISLQSL](#) (for pasting into other applications).
- Get peptides by:** ☒ Residue Number ☐ Increasing Mass ☐ Decreasing Mass
- Query:**

Start	End	Observed	Mr (exp1)	Mr (calc)	Delta	Score	Peptide
21	44	1769.1000	1768.4927	1768.9855	0.0072	9	QVNNNNR

MASCOT : Search Parameters

© 2007-2014 Matrix Science



Because glu-fib in SwissProt is not a tryptic peptide. The N-terminus is created by a post-translational cleavage after serine. If you now go back to the search form and select semi-trypsin or enzyme type none, you'll get the match.

Modifications

PMF ✓ SQ ✓ MS/MS ✓

Fixed modifications

Carbamidomethyl (C)

>

<

Acetyl (K)

Acetyl (N-term)

Acetyl (Protein N-term)

Amidated (C-term)

Amidated (Protein C-term)

Ammonia-loss (N-term C)

Biotin (K)

Biotin (N-term)

Carbamyl (K)

Carbamyl (N-term)

Carboxymethyl (C)

Variable modifications

Oxidation (M)

>

<

Display all modifications ☐

- Get details of current modifications, download updates, and define new entries at <http://www.unimod.org>
- User definable with an in-house Mascot installation

MASCOT : Search Parameters

© 2007-2014 Matrix Science

MATRIX SCIENCE

This screen shot shows how modifications are displayed in the search form in Mascot 2.3 and later. If you are using an earlier version, there are just two list boxes, one for fixed modifications and one for variable. In the current arrangement, you move modifications from the single list on the right to and from the lists on the left. This makes it easier to see at a glance what has been selected for the search. If the checkbox labelled 'Display all modifications' is clear, as shown here, you get a relatively short list of the most common modifications. If you check the box, a much longer list is available. You can keep your list of modifications up-to-date by downloading the latest information from Unimod. If you have a modification which you don't want to share with others, you can add it to the local configuration file. We'll describe how to go about doing this in detail in the Mascot Server Administration talk.

Modifications

PMF✓ SQ✓ MS/MS✓

Modifications

- Fixed / static modifications cost nothing
- Variable / differential modifications are very expensive
- Use minimum variable modifications, especially for PMF

Maybe oxidation of M

Maybe alkylation of C

MASCOT : Search Parameters

© 2007-2014 Matrix Science



Modifications in database searching are handled in two ways. First, there are the fixed or static or quantitative modifications. An example would be the efficient alkylation of cysteine. Since all cysteines are modified, this is effectively just a change in the mass of cysteine. It carries no penalty in terms of search speed or specificity.

In contrast, most post-translational modifications do not apply to all instances of a residue. For example, phosphorylation might affect just one serine in a peptide containing many serines. These variable or differential or non-quantitative modifications are expensive in the sense that they increase the time taken for a search and reduce its specificity. This is because the software has to permute out all the possible arrangements of modified and unmodified residues that fit to the peptide molecular mass. As more and more modifications are considered, the number of combinations and permutations increases geometrically. A so-called combinatorial explosion.

Hence, it is very important to be as sparing as possible with variable modifications. Especially in a peptide mass fingerprint, where the increase in the number of calculated peptides quickly makes it impossible to find a statistically significant match.

Quantitation

PMF✗ SQ✓ MS/MS✓

Quantitation

•More later ...

MASCOT : Search Parameters

© 2007-2014 Matrix Science




Quantitation is the subject of a separate presentation.

Protein mass

PMF✓
SQ✗
MS/MS✗

Protein mass kDa

- Applied as sliding window because there is no guarantee that the database entry represents the processed protein
- Slows down the search
- Never useful for MS/MS search. Only useful for Peptide Mass Fingerprint when
 - Analyte is small fragment of very large entry
 - Low complexity entry.


MASCOT : Search Parameters
© 2007-2014 Matrix Science


The protein mass is the mass of the intact protein in kDa applied as a sliding window. That is, the mass of the contiguous stretch of sequence which contains all of the matched peptide mass values. This will generally be less than the mass of the entire sequence entry. Consequently, if you specify a value for the protein mass, this acts only as a ceiling. Not only will you see smaller proteins on the hit list, you will also see larger ones, but all of the reported matches will be within a stretch of sequence less than or equal to the specified mass.

If this field is left blank, there is no restriction on protein mass

Specifying a protein mass will slow down the search a little.

It's hard to find examples where this parameter is useful. We include it mainly because many people requested it. It could give a better score if the analyte was small fragment of very large entry, or a low complexity protein. But, you can't know this in advance, so our general recommendation is to leave the protein mass open



This is the error window on experimental peptide mass values, not the error window for MS/MS fragment ion mass values, which is set using the MS/MS tol. \pm parameter.

Units can be selected from: percentage, milli-mass units, parts per million, or Daltons.

Specifying too tight a tolerance is a very common reason for failing to get a match.

Making an estimate of the mass accuracy doesn't have to be a guessing game. Protein View includes a graph of the mass errors for intact peptides. Just search a strong standard and look at the error graph. You'll normally see some kind of trend. Add on a safety margin and this is your error estimate. If you see something that looks like this, a mass tolerance of ± 0.5 Da is about right. It gives some safety margin. Remember that there will always be the odd outlier, like the data point at the lower left. It is the general trend and distribution of the majority of the data points that is important.

For a peptide mass fingerprint, the score depends on the peptide tolerance. In an MS/MS search, this parameter has no effect on the ions score. However, it does affect the search time. The larger the tolerance, the longer the search will take.

Peptide tolerance ¹³C

PMF ✗

SQ ✗

MS/MS ✓

13C

0

Sometimes, peak detection chooses the 13C peak

The normal test for a precursor match is:
 $TOL > absolute(exp - calc)$

If this field is set to 1, the test will also succeed for
 $TOL > absolute(exp - calc - 1)$

If this field is set to 2, the test will succeed for the above two conditions, plus:
 $TOL > absolute(exp - calc - 2)$

MASCOT : Search Parameters

© 2007-2014 Matrix Science

MATRIX SCIENCE

Sometimes, peak detection chooses the 13C peak rather than the 12C. In extreme cases, it may pick the 13C2 peak. The normal test for a precursor match is:

$TOL > absolute(exp - calc)$


Assuming the mass values and tolerance are in Da, if this field is set to 1, the test will also succeed for

$TOL > absolute(exp - calc - 1)$

If this field is set to 2, the test will succeed for the above two conditions, plus:

$TOL > absolute(exp - calc - 2)$

This means that you can use a tight mass tolerance and still get a match to a 13C peak. If you are using a very high accuracy instrument, note that the precise shifts are the carbon isotope spacings of 1.00335 and 2.00670, rather than 1 and 2.




This is the error window on MS/MS fragment mass values.

Units can be milli-mass units, Daltons, or ppm (Mascot 2.5 and later).

Specifying too tight or too loose a mass tolerance will reduce the ions score. Peptide View includes a graph of the mass errors for fragment ions.

Here, the mass tolerance is much too high. A more appropriate tolerance might be ± 0.3 . Having a tolerance which is much too high can sometimes lead to artefacts and false positives



Mass type specifies whether the experimental mass values are average or monoisotopic. Monoisotopic mass is the mass of the peptide where all atoms are the most abundant natural isotopes of their elements, e.g. Carbon 12, Nitrogen 14, Hydrogen 1, etc. In most cases, this is the first peak of the natural isotope distribution. Average mass is the chemical mass, which is the centre of gravity of the isotope distribution.

In Mascot, you cannot mix the two, and have (say) average precursors and monoisotopic fragments.

Most modern instruments produce monoisotopic mass values. You will only have an average mass if the entire isotope distribution has been centroided into a single peak, which usually implies very low resolution. If you get this setting wrong, the mass errors will be very large and show a strong trend, because the difference between an average and a monoisotopic mass for peptides and proteins is approximately 0.06%.

Charge

PMF✓
SQ✓
MS/MS✓


Mass values
☒ MH⁺
☐ M_v
☐ M-H⁻

Peptide charge
1+

- 1+ means MH⁺, 1- means M-H⁻, etc.
- For MS/MS, this setting is a default, which is rarely used.

MASCOT : Search Parameters

© 2007-2014 Matrix Science



These fields are used to specify the peptide charge state. The radio buttons are from the peptide mass fingerprint form. The drop down list is used on the sequence query and MS/MS forms.

The notation "1+", "2+", etc. is used to save space and because some HTML form fields do not support the use of superscripts and subscripts. "1+" always means MH⁺, "1-" always means M-H⁻, etc.

For MALDI-PSD, the precursor peptides will generally be MH⁺, so the charge state should be set to "1+"

For an MS/MS search, the value specified here is a default. Most peak lists always specify a charge state, so default is never used.

Data (PMF)

PMF✓ SQ✓ MS/MS✗

Data file

Query
NB Contents of this field are ignored if a data file is specified.

Browse...

- Mass [intensity] [additional text]
- Applied Biosystems Data Explorer (.pkm)
- Bruker Analysis AutoXecute Data Report
- Bruker XML
- mzData (1.05)
- mzML

MASCOT : Search Parameters

© 2007-2014 Matrix Science

The contents of the query window on the peptide mass fingerprint form are only used when no data file has been specified.

The data format for a peptide mass fingerprint is auto detected. It can be a simple list of mass values, one per line. If a second values is present, it is assumed to be intensity. Any further values on the same line are ignored

Mascot also supports other peak list formats, as listed.

mzData is the standard interchange format sponsored by the HUPO Proteomics Standards Initiative working group

Data (MS/MS)

PMF ✖
SQ ✖
MS/MS ✔

Data format

Mascot generic ▼


Precursor

m/z

- Mascot Generic Format (.MGF)
- Finnigan (.ASC)
- Sequest (.DTA)
- PerSeptive (.PKS)
- Micromass (.PKL)
- Sciex API III
- Bruker (.XML)
- mzData (.XML)
- mzML (.mzML)

MASCOT : Search Parameters

© 2007-2014 Matrix Science



Data for MS/MS ion searches must be supplied as an ASCII file in one of these supported formats. The format cannot be auto-detected, and must be specified using the drop down list.

Certain data file formats, SCIEX API III, PerSeptive (.PKS), and Bruker (.XML), do not include m/z information for the precursor peptide. For these formats only, the Precursor field is used to specify the m/z value of the parent peptide.

A data file may include embedded search parameters. Most embedded parameters can only appear once, at the head of the data file. In a Mascot generic format file, a few parameters can appear within an MS/MS dataset. See the Data File Format help page for further details

If there is a conflict between the values of the embedded parameters and values entered into search form fields, the embedded parameters always take precedence. The search form fields are essentially defaults for values missing from the data file.

Data URL

PMF ✗ SQ ✗ MS/MS ✓

☐ Data file

Choose file

No file chosen

Data input

☒ Data URL (http or ftp)

Data format


Mascot generic

Precursor

m/z

MASCOT : Search Parameters

© 2007-2014 Matrix Science

MATRIX SCIENCE


In Mascot 2.5 and later, if security is enabled, it is possible to specify a URL to the peak list. This means that the peak list file doesn't have to be downloaded to the client PC then uploaded to the Mascot server, which is useful for very large peak lists or when the client network connection is slow.

Instrument

PMF✗ SQ✓ MS/MS✓

Instrument ESI-QUAD-TOF

- Click on the help link to see which ions series are used



MASCOT : Search Parameters

© 2007-2014 Matrix Science

MATRIX
SCIENCE

For an MS/MS Ions Search, choose the description which best matches the type of instrument used to acquire the data. This setting determines which fragment ion series will be used for scoring, according to the following table.

Instrument

PMF ✗ SQ ✓ MS/MS ✓

MASCOT : Search Parameters

© 2007-2014 Matrix Science



"Default" corresponds to the configuration used in Mascot version 1.7 and earlier.

Many of the instruments are very similar.

You can modify instrument settings or create new ones using the configuration editor. In this screenshot, the right hand column is an experiment to see how the addition of w ions affects ETD matching

Error Tolerant

PMF ✗ SQ ✗ MS/MS ✓

Error tolerant ☐

The screenshot shows the Mascot search results page. At the top, there is a search bar with the text "Error tolerant" and a checkbox. Below the search bar, there is a table of results. The table has columns for "Query", "Observed", "Theoretical", "Delta Mass", "Score", "Expect", "Rank", and "Peptide". The results are sorted by score, with the highest score at the top. The table shows a list of peptides with their observed and theoretical masses, scores, and sequences. The sequences are shown in a monospace font, with some characters highlighted in red. The table is scrollable, and the bottom of the page shows the "MASCOT : Search Parameters" section.

Query	Observed	Theoretical	Delta Mass	Score	Expect	Rank	Peptide
1	1011.5000	1011.5000	0.0000	470	0.00	1	R.1000P1000.1
2	1011.5000	1011.5000	0.0000	470	0.00	2	R.1000P1000.1
3	1011.5000	1011.5000	0.0000	470	0.00	3	R.1000P1000.1
4	1011.5000	1011.5000	0.0000	470	0.00	4	R.1000P1000.1
5	1011.5000	1011.5000	0.0000	470	0.00	5	R.1000P1000.1
6	1011.5000	1011.5000	0.0000	470	0.00	6	R.1000P1000.1
7	1011.5000	1011.5000	0.0000	470	0.00	7	R.1000P1000.1
8	1011.5000	1011.5000	0.0000	470	0.00	8	R.1000P1000.1
9	1011.5000	1011.5000	0.0000	470	0.00	9	R.1000P1000.1
10	1011.5000	1011.5000	0.0000	470	0.00	10	R.1000P1000.1

MASCOT : Search Parameters

© 2007-2014 Matrix Science



If you have MS/MS data, and are interested in finding post-translational modifications, you can perform an error tolerant search by checking this box on the search form. This is a much more efficient way to discover unusual modifications, as well as non-specific peptides and sequence variants. More about this in a later presentation.

Decoy
PMF✓ SQ✓ MS/MS✓

☐ Decoy

The screenshot shows the 'Molecular & Cellular Proteomics' website. A table is highlighted with a red border, showing peptide matches above identity and homology thresholds. The table has columns for 'Sprot', 'Decoy', and 'False discovery rate'.

	Sprot	Decoy	False discovery rate
Peptide matches above identity threshold	3290	8	0.24 %
Peptide matches above homology or identity threshold	6037	224	3.71 %

MASCOT : Search Parameters
© 2007-2014 Matrix Science
 MATRIX SCIENCE

The decoy checkbox enables you to validate the false discovery rate according to the approach recommended in the Molecular & Cellular Proteomics Guidelines for Publication: “For large scale experiments, provide the results of any additional statistical analyses that indicate or establish a measure of identification certainty, or allow a determination of the false-positive rate, e.g., the results of randomized database searches or other computational approaches”

Report

PMF✓ SQ✓ MS/MS✓

Report top | AUTO hits

Report top should normally be set to auto.

MASCOT : Search Parameters© 2007-2014 Matrix ScienceMATRIX SCIENCE

REPORT determines the *maximum* number of hits displayed in a search results report. Choose AUTO to display only protein hits with significant scores. In a protein summary report, one additional hit is reported after the cutoff at the significant score. This is to ensure that the report shows the highest scoring hit, even though it is not significant.

Setting defaults

PMF✓ SQ✓ MS/MS✓

MASCOT : Search Parameters

© 2007-2014 Matrix Science

MATRIX SCIENCE

You can choose your own defaults for the search forms. Look for the link at the bottom of the search form selection page

When you save the defaults, they are saved as a browser cookie. If you go to a different PC, or switch to a different browser, you'll need to repeat this step

Final Tip

DANGER!

- Iteratively adjusting search parameters to get a better score can give misleading results
- Beware of
 - Narrowing the taxonomy
 - Reducing mass tolerances
 - Removing modifications
 - Selecting spectra or mass values

Set search parameters using standard samples

A final word of advice: It is easy to distort the search results without realising.

Basically, it is risky to adjust the search parameters interactively to get a better score for an unknown.

For example, you search the complete database and don't get a significant match. However, a very interesting looking protein is near the top of the list, surrounded by some others that are clearly wrong. You change the taxonomy filter so as to exclude the "wrong" proteins. Sorry, but this is cheating.

Search parameters should be set using standards. Broadening the search if you get a negative result is usually OK, but not narrowing the search.